

Monitoring of baker's yeast fermentations

by PD Dr.-Ing. Lars M. Blank, Laboratory of Chemical Biotechnology, Technical University Dortmund

● ● ● ● The open question we addressed with the new setup from BlueSens (CO₂ and ethanol sensor) originated from our previous finding (Blank and Sauer, 2004) that under aerobic and glucose excess conditions ethanol production and the rate of TCA cycle operation were dependent on the glucose uptake rate. As ethanol generally cannot be quantified in shake flasks, the finding relied only on indirect observations from ¹³C-tracer metabolic flux analyses. Here we aimed to directly quantify the TCA cycle flux by closing the carbon balance using the BlueSens sensors for quantification of the volatile fermentation products ethanol and CO₂. As can be seen in figure 1, the new setup delivered fer-


i

tu technische universität dortmund

PD Dr.-Ing. Lars M. Blank, Laboratory of Chemical Biotechnology, Technical University Dortmund. The group Systems Biotechnology characterizes, designs and constructs metabolic networks.

www.bci.tu-dortmund.de/bt

mentation data of very high quality (lines represent a simultaneous fit of the experimental data using an exponential growth model). As contribution to the scientific discussion, a strong negative correlation between glucose uptake rate and the rate of TCA cycle operation could be communicated (Heyland et al., 2009). The BlueSens setup was invaluable for the here presented quantitative physiology project with baker's yeast. Since then, numerous co-workers used the setup and experienced a tremendous increase in data amount and more importantly in quality.



Shake flasks equipped with CO₂, O₂ and ethanol sensor in a waterbath shaker




Figure 1. Fermentation course of *S. cerevisiae* during respiro-fermentative growth. (a) CO₂ and gaseous ethanol concentrations were monitored in the gas phase using infrared sensors. (b) Biomass plotted vs. CO₂ and gaseous ethanol concentrations. (c) Concentrations of glucose, ethanol, glycerol, and acetate were quantified by UV-RI-HPLC. (d) Biomass plotted vs. concentrations of glucose, ethanol, glycerol and acetate. Lines represent a best fit of all experimental data using an exponential growth model or by linear fit implemented in the Sigma Plot statistic module during exponential growth until 10 h. Linear fitting for gaseous CO₂ and Ethanol was only conducted until 9 h.

Literature

Blank, L. M. and U. Sauer, TCA cycle activity in *Saccharomyces cerevisiae* is a function of the environmentally determined specific growth and glucose uptake rate, *Microbiol.* 2004 150: 1085-1093

Heyland J., J. Fu, and L. M. Blank, Correlation between TCA cycle flux and glucose uptake rate during respiro-fermentative growth of *Saccharomyces cerevisiae*, *Microbiology*, 2009, 155: 3827-3837