



## Report from

P. Bakonyi, N. Nemestóthy, K. Bélafi-Bakó  
Research Institute on  
Bioengineering, Membrane  
Technology and Energetics  
University of Pannonia, Egyetem  
street 10, Veszprém, Hungary



## Application of BlueSens® H<sub>2</sub> sensor in a gas separation study

### Introduction

Biohydrogen is an emerging alternative energy carrier with high potential for sustainable development. It can be utilized in fuel cells or internal combustion engines. However, purified hydrogen should be provided for these end-use applications thus its enrichment – through getting rid of gaseous side-products especially CO<sub>2</sub> generated in the fermentative bioreaction – is very important.

A wide range of absorptive, adsorptive, cryogenic techniques and alternatives relying on membranes may be employed to concentrate biologically produced hydrogen. In recent decades as a result of advancements made in the field of membrane technology it is now accounted among the efficient, energy-saving and ecological-friendly opportunities. Membranes for gas separation are mature, commercially available nowadays and mostly fabricated from high-performance polymers such as polyimides. Membranes for H<sub>2</sub> purification, however, are featured with challenges. On one hand, membranes possess sufficiently high permeability and selectivity should be developed. On the other hand,

revealing the membranes real potential under practical separation conditions e.g. with mixed gases is still of interest but such studies are not frequently conducted despite their importance.

Therefore, in this work we wanted to test a lab-scale polyimide membrane module for fermentative hydrogen recovery applying binary H<sub>2</sub>/CO<sub>2</sub> gaseous mixtures.

Moreover, evaluating the influences of notable process variables, namely separation temperature and recovery factor was also aimed.

### Materials and Methods

#### Hydrogen purification experiments using BlueSens® H<sub>2</sub> sensor

In the current experiments a membrane module containing polyimide hollow-fibers (UBE Industries Ltd., NM-B01A) was investigated. This polymer is considered as a highly promising material for membrane gas separation including H<sub>2</sub> purification, as well.

The measurements aiming the concentration of hydrogen were carried out in a laboratory-scale device schematically presented in in Fig. 1.




Fig. 1 – Experimental membrane testing set-up  
1; gas cylinder containing binary H<sub>2</sub>/CO<sub>2</sub> mixture, 2; pressure controller, 3; feed valve, 4; feed pressure indicator, 5; module thermostat, 6; membrane module, 7; differential pressure indicator, 8; permeate pressure indicator, 9,10; permeate valves, 11,15; digital gasflow meter, 12,16; BlueSens® H<sub>2</sub> sensor, 13; permeate fraction, 14; retentate valve, 17; retentate fraction

During the experiments, the hydrogen concentrations of the feed, permeate and retentate streams were continuously followed by using on-line BlueSens® H<sub>2</sub> analyzers (BlueSens Gas Sensor GmbH). The equipment is able to measure hydrogen in the range of 0-100 vol% with a quick response time of 20 seconds. The additional characteristics can be accessed by visiting the manufacturer's homepage ([www.bluesens.com](http://www.bluesens.com)).

The feed gas pressures during the experiments were adjusted as low as approximately 3 bar since considering practicality, a high compression degree of biologically produced gases comprising hydrogen would consume a noticeable portion of its energy content and hence making the process economically unattractive. The permeate side of the module was kept under ambient conditions.

The influences of three process variables (temperature, feed hydrogen content and retentate/feed flow ratio referred to as recovery factor) were determined by experimental design approach using binary H<sub>2</sub>/CO<sub>2</sub> mixtures.

To estimate the mixed gas selectivities ( $a_{i/j}$ ) according to Eq. 1, only the data recorded at steady-state – when no change in the retentate and permeate flows and in their compositions occurred – were considered:

$$a_{i/j} = (x_i^P/x_j^P)/(y_i^F/y_j^F) \quad (1)$$

where  $x_i^P$  and  $x_j^P$  are the volumetric fractions of compound „i” and „j” in the permeate, meanwhile  $y_i^F$  and  $y_j^F$  are the their volumetric fractions in the feed mixture.

## Results and discussion

Statistical analysis (ANOVA) was performed to evaluate the experimental design (Table 1). Based on ANOVA results, feed H<sub>2</sub> concentration, temperature and recovery factor were all identified as significant process parameters ( $p < 0.05$ ).

The response plots displaying the impacts of the independent variables mentioned can be seen in Fig. 2-Fig. 4.

### Impact of gas composition

The composition of the feed gas is apparently a factor that determines the performance of the membrane module. However, the separation efficiencies of the membranes are more often than not sought only in single gas experiments. This ideal value can be notably different to those obtained under practical (mixed gas) separation circumstances. The theoretical and real selectivities would expectedly be similar only when the constituents of a gas mixture do not strongly interact with the membrane material. In the case of CO<sub>2</sub> and H<sub>2</sub> separation, however, the carbon-dioxide is easily dissolved in membrane materials such as polyimide. This phenomenon can remarkably affect the permeation of the other gases e.g. hydrogen and thus

the efficiency of the overall separation process. From the tentative results, we found that higher hydrogen concentration could be attained in the permeate fraction by feeding  $H_2/CO_2$  mixtures with higher initial hydrogen concentration. Based on the experiments it was concluded that the more amount of carbon-dioxide in the feed mixture represented a limiting factor to achieve better selectivities. This might be explained by the fact that carbon dioxide likely influences the migration of hydrogen across the membrane capillaries. As  $CO_2$  concentration has decreased, improved separation efficiency could be observed.

## ANOVA – Dependent Variable: Mixed gas selectivity;

**R2: 0.98**

Factor	SS	df	MS	F	p
Separation Temperature	0.007200	1	0.007200	11.2941	0.02828
Feed $H_2$ content	0.068450	1	0.068450	107.3725	0.00049
Recovery factor	0.051200	1	0.051200	80.3137	0.00086
Error	0.002550	4	0.000637		




Fig. 2 – Influence of separation temperature and feed composition on mixed gas selectivity




Fig. 3 – Influence of feed composition and recovery factor on mixed gas selectivity






Fig. 4 – Influence of separation temperature and recovery factor on mixed gas selectivity

#### Influence of recovery factor

Membranes, in general, can enrich the slowly permeating compounds in the concentrate stream, however, only when low recoveries are applied. In our investigation, almost pure (>99 vol%) CO<sub>2</sub> could be obtained in the retentate at certain, extremely low recoveries when the vast majority of the gas introduced to the membrane module was taken at the permeate side (Table 2).

CO <sub>2</sub> Permeate (vol%)	CO <sub>2</sub> Retentate (vol%)	Recovery value	Mixed Gas Selectivity
67.17	99.61	0.05	1.12
65.42	83.8	0.17	1.21
63.82	76.13	0.38	1.3
62.8	71.74	0.69	1.36

Table 2 – Results of separating H<sub>2</sub>/CO<sub>2</sub> gas mixture (Separation temperature: 37 °C; initial H<sub>2</sub> content: 30 vol%)

It was also observed that the separation of the faster permeating H<sub>2</sub> is determined by the membrane's selectivity, since H<sub>2</sub> concentration in the permeate was only slightly increased when recovery factor was increased. Hence, if hydrogen is to be concentrated in the permeate fraction, membranes reflecting higher

selectivities need to be developed or multi-step, cascade applications should be employed. Furthermore, the change in the recovery value also changes the amount of permeate which together with the total feed flow and the surface velocities should not be neglected to avoid concentration polarization. Basically, maintaining a high recovery value is favorable so as to get a more sufficient membrane performance. This might be the reason for the increment of selectivity with increased recoveries.


#### Separation temperature effect

It is quite obvious from the research outcomes that the mixed gas selectivity has increased at higher temperatures. It is understandable if one keeps in mind that increasing temperature has contradictory impacts on the solubility and diffusivity of the permeating gases. In general, the former decreases, whilst the latter increases with raising temperatures. Furthermore, hydrogen and carbon dioxide have distinct characteristics from diffusion and solubility points of views in rigid, glassy polymers such as polyimide. CO<sub>2</sub> usually expresses a relatively higher solubility in comparison to hydrogen which is in turn more rapid and possesses faster diffusion coefficient. These properties of the gases are in correlation with molecular size, their affinity to polymers and the developing interactions with the membrane material during permeation. Polyimide, as mentioned, is a glassy polymer which basically achieves selectivity mainly on diffusivity-difference basis. Since CO<sub>2</sub> is a larger molecule than hydrogen, the increment in its diffusivity is less pronounced than that of hydrogen at elevated temperatures. This was assumed to be the reason for the increasing separation efficiency of H<sub>2</sub>/CO<sub>2</sub> with increasing temperature.

#### Conclusions

In our research work we established that the polyimide membrane module had a potential to enrich hydrogen from gaseous mixtures. The results demonstrated that process variables such as gas composition, temperature and the recovery factor had remarkable

influence on hydrogen separation and could affect the achievable separation efficiency. Thus, these parameters are to be considered when designing a suitable enrichment process. It has been proven that higher  $H_2/CO_2$  mixed gas selectivity could be accomplished by increasing the feed hydrogen concentration as well as the temperature and the recovery value within the studied design boundaries.



Dipl. Eng Péter Bakonyi, University of Pannonia  
Research Institute on Bioengineering, Membrane  
Technology and Energetics

Ph.D. Nándor Nemestóthy, University of Pannonia  
Research Institute on Bioengineering, Membrane  
Technology and Energetics

Prof. Katalin Bélafi-Bakó, University of Pannonia  
Research Institute on Bioengineering, Membrane  
Technology and Energetics

The main tasks of the Research Institute on  
Bioengineering, Membrane Technology and Energetics  
are:

- > to teach and train students to get BSc In  
bioengineering
- > to introduce them to a wide range of  
biochemical and biotechnological processes
- > to introduce them to the world of membranes  
and their applications

Application possibilities of various membrane  
processes, improvements of bioprocesses and  
manufacture of renewable, “green” energy sources are  
studied at the institute involving the PhD students of  
the Doctoral School of Chemical Engineering and  
Material Sciences (University of Pannonia).