English

Deutsch

# Operating Manual

## BCP-CH4

**BlueSens**

• • • • • • • •

## Contents

1	ABOUT THIS DOCUMENT .....	3
1.1	Function.....	3
1.2	Target group.....	3
1.3	Symbols used.....	3
2	FOR YOUR SAFETY.....	4
2.1	General information.....	4
2.2	Authorized personnel .....	4
2.3	Proper use.....	4
2.4	Misuse warning .....	4
2.5	General safety information .....	4
2.6	CE conformity.....	5
3	PRODUCT DESCRIPTION .....	5
3.1	One-piece construction of the BCP-CH4 .....	5
3.2	Measuring principle .....	5
4	INSTALLATION.....	6
4.1	General instructions .....	6
4.2	Mechanical connection.....	6
4.2.1	Installation on pipes .....	7
4.2.2	Installation on a Tri-Clamp SMS38 connection.....	7
4.2.3	Installation on a POM flow adapter .....	8
4.2.4	Installation on a stainless steel flow adapter.....	8
4.2.1	Installation on pipes .....	9
4.2.2	Installation outside .....	9
4.2.3	Ball valve .....	9
4.3	Electrical connection .....	10
4.3.1	General information.....	10
4.3.2	Version 4 – 20 mA in PA6 housing .....	11
4.3.3	RS232 serial version in PA6 housing .....	12
4.3.4	Version 4 – 20 mA in aluminum housing .....	13
4.3.5	RS232 serial version in aluminum housing .....	14
4.3.6	Connection via BACCom12 .....	15
4.4	Minimization of dilution effects .....	16
4.5	Cross-sensitivity to CO <sub>2</sub> and hydrocarbons.....	17
5	MAINTENANCE .....	18
5.1	1-point calibration .....	18
5.2	Recalibration .....	18
5.3	Filter change – coarse filter .....	19
5.3.1	Removing the filter cover .....	19
5.3.2	Changing the filter and seals (Z-XX-00052).....	19
6	EXTERNAL PRESSURE SENSOR P2518G AND P2518NPT .....	20
7	APPENDIX .....	22
7.1	Calibration table .....	22
7.2	Technical data .....	23

## 1 About this document

### 1.1 Function

This operating manual provides you with all of the necessary information for quick start-up and safe operation of the **BCP-CH4**. Please read the operating manual before starting operation.

### 1.2 Target group

This operating manual is intended for use by trained specialist personnel. The contents of this manual must be made available to personnel and followed by them.

### 1.3 Symbols used

#### Danger!

This symbol indicates a situation that is possibly dangerous. Failure to observe the safety instructions can result in personal injury.

#### Caution!

The symbol indicates the possibility of damage to property.

#### Note!

This symbol indicates helpful additional information.

#### ● List

This symbol indicates a list. The order of the items is not significant.

#### 1 Action sequence

Numbers indicate steps to be performed in sequence.

## 2 For your safety

### 2.1 General information

The **BCP-CH4** was inspected in our plant and was ready for operation when it left.

Before installing and starting up the device, please read this operating manual carefully. The operating manual contains safety instructions that must be observed to ensure safe operation.

The device must never be operated in conditions that do not comply with the specifications on the type plate.

Maintenance and servicing may only be performed by specially trained personnel who are familiar with the hazards inherent to the work as well as the guarantee terms.

### 2.2 Authorized personnel

All of the actions described in this operating manual may only be performed by trained specialist personnel who have been authorized by the plant operator. Work on the device other than that described in this manual may only be performed by personnel of the BlueSens gas sensor GmbH Company for safety reasons and to ensure compliance with the terms of the guarantee.

### 2.3 Proper use

The **BCP-CH4** is a gas sensor for measuring methane concentrations in the specified concentration area and under the conditions described in the technical data. It is used to monitor metabolism in biological processes such as fermentation. The BCP-CH4 sensor may only be used in well ventilated rooms.

#### Danger!

The sensor does not have an ATEX certificate and may therefore only be used in well ventilated rooms.

### 2.4 Misuse warning

The **BCP-CH4** may not be used as a safety component for monitoring gasses in systems or as a gas warning device. It may also not be used in areas subject to explosion hazards.

### 2.5 General safety information

If the device is mishandled or not used for its intended purpose, application-specific dangers may arise.

#### Danger!

If the device is incorrectly installed or set, there is a danger of explosions and poisoning.

After installation, check all connections for leaks.

## 2.6 CE conformity

The **BCP-CH4** conforms to the EMC Directive (89/336/EEC, 92/31/EEC and 93/68/EEC) when applying the harmonized standards **EN50081-1, EN61000**.

The low-voltage directive (72/23/EEC und 93/68/EEC) is not applicable as no voltage greater than 24 V is used.

## 3 Product description

### 3.1 One-piece construction of the BCP-CH4

The one-piece construction (fig. 1) means that the measuring adapter cannot be separated from the sensor head. It is the standard set-up for all BCP-sensors. The BCP-CH4 is designed for a particular mechanical connection that can only be altered subsequently at the plant for a certain fee.


*Fig. 1: One-piece construction of the BCP-CH4*

### 3.2 Measuring principle

The sensor contains the IR light source, the detector and the evaluation electronics (fig. 2).

The infrared light beam is reflected by the gas-filled measuring adapter and the light weakened by the analyte gas is measured by the detector. The light-permeable sapphire disc prevents the sample atmosphere from escaping and contaminating the optical components. The sensor head heats the measuring adapter so that no moisture can condensate. The heating time lasts about 1 hour both during initial start-up and after each time the device is disconnected from the power supply. The sensor does not output any measurement values during the heating-up period.

*Fig.2: Measuring principle*

A	Sensor head
B	Detector
C	Light source
D	Seal
E	Sapphire window
F	Light beam
G	Sterile filter

## 4 Installation

### 4.1 General instructions

The **BCP-CH4** is protected by packaging on its way to its application location. This secures it against the usual transport strains. However, before installation, check whether the device has been damaged due to improper transport or improper storage. If the device is damaged in any way, operation without hazards is not possible and the device may not be installed and taken into operation.

Check whether the enclosed materials such as seals and screw-caps are suitable for your process conditions (pressure, temperature, etc.).

The installation should only be performed under supervision by a specialist and in compliance with all applicable work safety rules.

### 4.2 Mechanical connection

#### Caution!

The sterile filter is not intended for repelling fluids. Never install the sensor such that fluid can run into the measuring adapter.

If water has penetrated the measuring adapter, allow it to dry out for at least 12 hours at max. 80°C in a drying cabinet or on a hot plate.

Protect the measuring adapter from penetration by liquids.

After installation, check that the pipe connection is gas-tight.

#### 4.2.1 Installation on pipes

The connection to a pipe is made with a 1 ¼" nozzle with an external thread:

1. Place the sealing ring (30 x 4 mm viton, item no. Z-OR-00003) on the nozzle (fig. 3).
2. Place on the sensor (fig. 4).
3. Connect the nozzle and the sensor with the screw cap so that the connection is gas-tight (fig. 5).

##### Note!

**Only use the supplied screw caps. Do not use metal screw caps as they result in thermal contact between the measuring adapter and the pipe and thus violate the technical specifications.**

Fig. 3

Fig. 4

Fig. 5

#### 4.2.2 Installation on a Tri-Clamp SMS38 connection

:

1. Place the sealing ring (item no. Z-OR-00013) on the nozzle (fig. 6).
2. Place on the sensor (fig. 7).
3. Fix the sensor with the Tri-Clamp on the nozzle (fig. 8).

Fig. 6

Fig. 7

Fig. 8

#### 4.2.3 Installation on a POM flow adapter

To install the sensor on a POM flow adapter:

1. Place the sealing ring (item no. Z-OR-00004) on the nozzle of the flow adapter (fig. 9).
2. Place on the sensor (fig. 10).
3. Connect the flow adapter and the sensor with the screw cap so that the connection is gas-tight (fig. 11).

Fig. 9

Fig. 10

Fig. 11

#### 4.2.4 Installation on a stainless steel flow adapter

To install the sensor on a stainless steel flow adapter:

1. Place the sealing ring (Z-OR-00004) on the stainless steel connection piece (fig. 12).
2. Place on the sensor and put the screws in place (fig. 13).
3. Fasten the 4 screws (item no. Z-XX-00007) so that the connection is gas-tight (fig. 14).

Fig. 12

Fig. 13

Fig. 14

#### 4.2.5 Installation on pipes

The ball valve can be connected to the 1 ¼" pipe via a screwed connection.

1. Remove the protective plastic cap from the valve (fig.15).
2. Place the valve on the screw thread of the pipe (fig.16).
3. Use Locktite™, hemp-sealing, Teflon-band or other thread sealant to achieve a gas-tight connection.
4. Tighten the screwed connection.

#### Note!

**After installation, check for gas-tightness.**

#### 4.2.1 Installation outside

The **BCP-CH4** complies with the Ingress Protection (IP) class 65. The sensor can be used outside for a long time. Under extreme climatic conditions like permanent rain and storm or cold temperatures we strongly recommend to protect the sensor with an additional protective housing against influences from the weather. This housing could consist of a cubic steel housing with two open sides (Fig 17).

#### 4.2.2 Ball valve

With the ball valve it is possible to disconnect the sensor from the gas stream for maintenance. During the normal measurement the ball valve is in the "On"-position (fig.18 so that the analyte gas and the pressure can be detected by the BCP-CH4).

For maintenance, like the 1-point-calibration (see chapter 5.1) the sensor can be cut off the gas stream. Pull the handle of the valve on the position "off" (fig. 19) to disconnect the sensor from the gas stream. The process can continue without measurement in the meantime when indicated.

Fig. 15: Removal of the plastic cap


Fig. 16: Valve on screw thread

Fig. 17: additional protection


Fig. 18 valve open

Fig. 19 valve closed

## 4.3 Electrical connection

### 4.3.1 General information

#### Caution!

Read the installation instructions carefully to prevent damage to the device.

Proceed step-by-step.

Only use the original plugs, cables and power adapters.

Never connect or disconnect plugs when the device is connected to the power supply.

The device does not have an on/off switch; it starts operation as soon as it is connected to the power supply.

Improper operation can result in damage to the device.

#### Note!

New wide temperature sensor in aluminum housing (WT sensors) need a power supply with a voltage of 24V. They will not work with 12V. Please check our individual date sheet.

### 4.3.2 Version 4 – 20 mA in PA6 housing

To connect the measuring device to the connection cable of the sensor head in the PA6 housing (fig. 20), use the supplied socket and strain relief (fig. 21).

#### Note!

**The numbering of the pins and their assignment refer to the socket when seen from behind (fig. 22).**

**Remove the insulation of the cables a little as possible to avoid short circuits in the plug housing.**

Fig. 20

Fig. 21

PIN 1	V+ = 12 V
PIN 2	GND
PIN 3	RS232_TXD
PIN 4	RS232_RXD
PIN 5	1-point calibration
PIN 6	4–20mA, RL < 250 Ohm
PIN 7	For internal use only
PIN 8	GND



Fig. 22: Plug assignment

1. Connect the 12 V DC power supply to pin 1 of the socket.
2. Connect GND (ground) to pin 2.
3. Connect the measuring device to pin 6 ( $R_L < 250 \text{ Ohm}$ ) and pin 8 GND (ground).
4. Plug the sensor cable into the socket.

After around 1 hour of heating-up time, the sensor still requires adjusting. During the heating-up time, the sensor displays 2.3 mA. To make the adjustment, expose the sensor for approx. 30 minutes (depending on specification – see datasheet) to normal fresh air (0.0 vol. % CH4) or N<sub>2</sub>.

5. Afterwards, connect pin 5 to pin 8 (GND) for 5 seconds.
6. Screw on the strain relief. The sensor has been adjusted.

#### 4.3.3 RS232 serial version in PA6 housing

1. Connect the sensor to the power supply with the cable supplied.
2. Connect the sensor to a computer using the serial cable.

After around 1 hour of heating-up time, the sensor still requires adjusting. To do this, expose the sensor for approx. 30 minutes (depending on specification – see datasheet) to normal fresh air (0.0 vol. % CH<sub>4</sub>) or N<sub>2</sub>.

The adjustment itself is performed with the **BACVisSingle** software (see **BACVisSingle** operating manual).

Start the **BACVisSingle** software. You will find further relevant information in the corresponding operating manual.

Fig. 23

#### 4.3.4 Version 4 – 20 mA in aluminum housing

To connect the measuring device to the connection cable of the sensor head in the aluminum housing (fig. 24), use the supplied socket and strain relief (fig. 25).

##### Note!

**The numbering of the pins and their assignment refer to the socket when seen from behind (fig. 26).**

**Remove the insulation of the cables a little as possible to avoid short circuits in the plug housing.**


Fig. 24

Fig. 25

PIN 1	V+ see voltage on your individual data sheet
PIN 2	GND
PIN 3	RS232_TXD
PIN 4	RS232_RXD
PIN 5	1-point calibration
PIN 6	4–20 mA, $R_L < 250 \text{ Ohm}$
PIN 7	For internal use only
PIN 8	GND



Fig. 26: Plug assignment

1. Connect the DC power supply to pin 1 of the socket. See voltage on your individual data sheet.

2. Connect GND (ground) to pin 2.

3. Connect the measuring device to pin 6 ( $R_L < 250 \text{ Ohm}$ ) and pin 8 GND (ground).

4. Plug the sensor cable into the socket.

After around 1 hour of heating-up time, the sensor still requires adjusting. During the heating-up time, the sensor displays 2.3 mA. To make the adjustment, expose the sensor for approx. 30 minutes (depending on specification – see datasheet) to normal fresh air (0.0 vol. % CH4) or N<sub>2</sub>.

5. Afterwards, connect pin 5 to pin 8 (GND) for 5 seconds.

6. Screw on the strain relief. The sensor has been adjusted.

#### 4.3.5 RS232 serial version in aluminum housing

1. Connect the sensor to the power supply with the cable supplied. See voltage on your individual data sheet.
2. Connect the sensor to a computer using the serial cable.

After around 1 hour of heating-up time, the sensor still requires adjusting. To do this, expose the sensor for approx. 30 minutes (depending on specification – see datasheet) to normal fresh air (0.0 vol. % CH<sub>4</sub>) or N<sub>2</sub>.

The adjustment itself is performed with the **BACVisSingle** software (see **BACVisSingle** operating manual).

Start the **BACVisSingle** software. You will find all of the further relevant information in the corresponding operating manual.

Fig. 27

<b>1 = +12 or + 24 V, please see you data sheet</b>
<b>2 = 0 V</b>
<b>3 = RS232_RxD</b>
<b>4 = RS232_TxD</b>
<b>5 = RS232_GND = PE</b>
PE = ground

Fig. 28: Plug assignment

#### 4.3.6 Connection via BACCom12

The **BACCom12** connection box is an electronic multiplexer with an integrated pressure sensor. It facilitates the connection of up to 12 sensor heads.

Communication with a PC can be switched between RS232 or Ethernet.

The individual connections are explained in the following table:

	Designation	Description
A	RJ45	RJ45 socket for connecting the sensors
B	LED	Operating display when a voltage is present
C	Sub D 9 pin	Data transmission to the PC
D	Switch	Switches between RS232 and Ethernet
F	RJ45	Ethernet connection
G	Power socket	12 V 3.75 A, only use the supplied power adapter
H	Box reset	Resets the box; does not effect the sensors
K	M8 4 pin socket	4-pin connection sockets A–D for additional boxes

#### Caution!

To prevent damage to the device, only use the supplied power adapter and the supplied cable.

Never disconnect or connect the connection plugs on the sensor heads when the **BACCom12** is switched on.

1. Connect all sensor heads with the **BACCom12**.
2. Connect the supplied power adapter to the power socket **G**.
3. Plug the power plug of the power adapter into the socket.

After a heating-up time of approx. 1 hour, the measuring system is ready for operation.

Fig. 29: Front of the **BACCom12**

Fig. 30: Connections on the **BACCom12**

4. Connect the **BACCom12** to the PC or network via the Ethernet port **E**, or connect the **BACCom12** via the RS232 output **C** with the **supplied** cable to the serial interface of the computer.
5. Select the corresponding interface with the switch **D**.

After around 1 hour of heating-up time, the sensor still requires adjusting. To do this, expose the sensor for approx. 30 minutes (depending on specification – see datasheet) to normal fresh air (0.0 vol. % CH4) or N<sub>2</sub>.

Adjustment of the sensors is performed with the **BACVis** software. Start the corresponding software. You will find all further information in the software instructions.

After initial commissioning, the measuring device can remain switched on constantly, meaning that the heating-up time is not required before every measuring.

#### 4.4 Minimization of dilution effects

The dilution effect due to moistening of dry gas can be compensated by dividing the displayed concentration value by the dilution factor.

To do this, you need to know the temperature directly at the output of the exhaust cooler or the temperature in the gas line close to the sensor. Use the lower of the two temperatures to determine the dilution factor from the table. A linear regression can be assumed between the individual values.

Temperature [°C]	Dilution factor
-10	0.9972
0	0.9940
5	0.9914
10	0.9879
15	0.9831
20	0.9769
25	0.9687
30	0.9580
35	0.9444
40	0.9270
45	0.9052
50	0.8779

#### 4.5 Cross-sensitivity to hydrocarbons and CO<sub>2</sub>

Due to the NDIR-measurement method of the BCP-CH4, every sensor has got a principle-related cross-sensitivity to other hydrocarbons. The CO<sub>2</sub>'s absorption band is very close to the reference wavelength of the detector. Therefore there is also a certain cross-sensitivity to CO<sub>2</sub>.

The value of deviation caused by the cross-sensitivity is dependent on the individual measurement specification of the BPC-CH4 and also on the gas applied. For the measurement of biogas the BCP-CH4 Biogas is calibrated in a special way to minimize the cross-sensitivity to CO<sub>2</sub>. For more details according the cross-sensitivity and your individual sensor please contact us directly.

## 5 Maintenance

We recommend sending the device to BlueSens for annual maintenance, checking and calibration of the sensors.

### 5.1 1-point calibration

Once monthly, or after each connection and disconnection of the sensor head and measuring adapter, the sensor head must be exposed for approx. 30 minutes (depending on specifications, see datasheet) to normal fresh air (0.0 vol. % CH<sub>4</sub>) or nitrogen (N<sub>2</sub>).

Afterwards, connect pin 5 to pin 8 on the connection cable for 5 seconds, or, if present on the sensor, press the blue button for 5 seconds (fig. 31).

For the serial version, the adjustment can be made using the **BACVis** software.

Fig. 31

### 5.2 Recalibration

The sensor should be sent back to the manufacturer or an authorized dealer for annual recalibration.

You can get further information for our annual maintenance service Blue4Care incl. extension of the warranty up to 6 years on

<http://www.bluesens.de/fileadmin/dl/Blue4Care.pdf>

## 5.3 Filter change – coarse filter

### 5.3.1 Removing the filter cover

Hold the sensor head and screw off the lower cover counterclockwise using the clamping ring tongs (fig. 32). If your filter-cover is already equipped with notches (fig. 33) you can use a big coin (like a 2-Euro-coin or a quarter) to remove the cover, too (fig. 34).

Fig. 32

### 5.3.2 Changing the filter and seals (Z-XX-00052)

1. Remove the filter from the recess (fig. 35).
2. Insert a new filter.
3. Check the seals for damage and replace as required.
4. Screw on the lower cover using the clamping ring tongs respectively using the coin.

Fig. 33


Fig. 34

Fig. 35

## 6 External pressure sensor P2518G and P2518NPT

The external pressure sensor P2518G/P2518NPT measures the pressure directly in the gas line. Inaccurate results by changing pressures in the process can be prevented by this device. The pressure sensor is separately available for the BCP-sensors in aluminium housing. An external pressure sensor can be ordered with a new sensor or retrofit to an existing sensor in the factory. Using the ordered pressure sensor is mandatory as the BCP gas sensor needs the pressure signal for producing proper measuring results. Wrong measuring results will occur if you do not use both devices together after the BCP gas sensor was aligned to have an external pressure sensor connected.

The pressure sensor (fig. 36:1) has got a 1/8" NPT male screw thread. With the screw thread adapter (fig. 36:2) the pressure sensor can be adapted to a G 1/8" male screw thread. With these screw threads the pressure sensor can be integrated directly into the gas line. Please use an appropriate thread sealant to provide a gas tight connection. Check the connection for tightness afterwards.

*Fig.36: pressure sensor and adapter*

*Fig.37: connection of pressure sensor to gas sensor*

The pressure sensor is connected to the corresponding gas sensor via a 4-pin M12 male connector (fig. 37:3). Plug the male connector into the female connector and tighten the thread connection (fig. 37:4). Each pressure sensor is matched individually with the corresponding gas sensor. Therefore the combination of gas sensor and pressure sensor belongs together and should not be interchanged with other components.

The pressure sensor was factory-calibrated for a specified pressure range. Please see the corresponding data sheet for specification.

## 7 Appendix

## 7.1 Calibration table

Complete calibration can only be conducted by BlueSens. Monthly 1-point calibration can be performed as described in chapter 5.1. The adjustment must also be made each time the sensor head and the measuring adapter are disconnected from each other.

Fill out the table below when this is performed.

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We hereby declare that our product, type:

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complies with the following relevant provisions:

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Low voltage guideline (72/23/EEC and 93/68/EEC) is not applicable as no voltage higher than 24 V is used.

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Applied harmonized standards, in particular:

**EN50081-1**  
**EN61000**

A handwritten signature in blue ink, appearing to read "Udo Schmale".

Dr. Holger Mueller, Dr. Udo Schmale  
BlueSens gas sensor GmbH  
Snirgelskamp 25  
D-45699 Herten, Germany  
Phone +492366 / 4995-500  
Fax +492366 / 4995-599  
[www.bluesens.com](http://www.bluesens.com)

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