English

Deutsch

Operating Manual

Multilingual

BCP-H2

BlueSens

• • • • • • • •

Contents

1	ABOUT THIS DOCUMENT	3
1.1	Function.....	3
1.2	Target group.....	3
1.3	Symbols used.....	3
2	FOR YOUR SAFETY.....	4
2.1	General information.....	4
2.2	Authorized personnel	4
2.3	Proper use.....	4
2.4	Misuse warning	4
2.5	General safety information	4
2.6	CE conformity.....	5
3	PRODUCT DESCRIPTION	5
3.1	One-piece construction of the BCP-H2.....	5
3.2	Measuring principle	5
4	INSTALLATION.....	6
4.1	General instructions	6
4.2	Mechanical connection.....	6
4.2.1	Installation on pipes	7
4.2.2	Installation on a Tri-Clamp SMS38 connection.....	7
4.2.3	Installation on a POM flow adapter	8
4.2.4	Installation on a stainless steel flow adapter.....	8
4.3	Electrical connection	9
4.3.1	General information.....	9
4.3.2	Version 4 – 20 mA in PA6 housing	10
4.3.3	RS232 serial version in PA6 housing.....	11
4.3.4	Version 4 – 20 mA in aluminum housing	12
4.3.5	RS232 serial version in aluminum housing	13
4.3.6	Connection via BACCom12	14
4.4	Cross-sensitivity to other gases	16
5	MAINTENANCE	17
5.1	1-point calibration.....	17
5.2	Recalibration	17
5.3	Filter change – coarse filter	18
5.3.1	Removing the filter cover	18
5.3.2	Changing the filter and seals (Z-XX-00052).....	18
6	APPENDIX	19
6.1	Calibration table	19
6.2	Technical data	20

1 About this document

1.1 Function

This operating manual provides you with all of the necessary information for quick start-up and safe operation of the **BCP-H2**. Please read the operating manual before starting operation.

1.2 Target group

This operating manual is intended for use by trained specialist personnel. The contents of this manual must be made available to personnel and followed by them.

1.3 Symbols used

Danger!

This symbol indicates a situation that is possibly dangerous. Failure to observe the safety instructions can result in personal injury.

Caution!

The symbol indicates the possibility of damage to property.

Note!

This symbol indicates helpful additional information.

1 Action sequence

Numbers indicate steps to be performed in sequence.

2 For your safety

2.1 General information

The **BCP-H2** was inspected in our plant and was ready for operation when it left.

Before installing and starting up the device, please read this operating manual carefully. The operating manual contains safety instructions that must be observed to ensure safe operation.

The device must never be operated in conditions that do not comply with the specifications on the type plate.

Maintenance and servicing may only be performed by specially trained personnel who are familiar with the hazards inherent to the work as well as the guarantee terms.

2.2 Authorized personnel

All of the actions described in this operating manual may only be performed by trained specialist personnel who have been authorized by the plant operator. Work on the device other than that described in this manual may only be performed by personnel of the BlueSens gas sensor GmbH Company for safety reasons and to ensure compliance with the terms of the guarantee.

2.3 Proper use

The **BCP-H2** is a gas sensor for measuring hydrogen concentrations in the specified concentration area of binary gas mixtures and under the conditions described in the technical data.

Danger!

The sensor does not have an ATEX certificate and may therefore only be used in well ventilated rooms.

2.4 Misuse warning

The **BCP-H2** may not be used as a safety component for monitoring gasses in systems or as a gas warning device. It may also not be used in areas subject to explosion hazards.

2.5 General safety information

If the device is mishandled or not used for its intended purpose, application-specific dangers may arise.

Danger!

**If the device is incorrectly installed or set, there is a danger of explosions and poisoning.
After installation, check all connections for leaks.**

2.6 CE conformity

The **BCP-H2** conforms to the EMC Directive (89/336/EEC, 92/31/EEC and 93/68/EEC) when applying the harmonized standards **EN50081-1, EN61000**.

The low-voltage directive (72/23/EEC und 93/68/EEC) is not applicable as no voltage greater than 24 V is used.

3 Product description

3.1 One-piece construction of the BCP-H2

The one-piece construction (fig. 1) means that the measuring adapter cannot be separated from the sensor head. It is the standard set-up for all BCP-sensors. The BCP-H2 is designed for a particular mechanical connection that can only be altered subsequently at the plant for a certain fee.

Fig. 1: One-piece construction of the BCP-H2

3.2 Measuring principle

The hydrogen sensor is based on a micro-mechanical silicon chip with a thin membrane. Two thin film resistors are integrated into the membrane and are used for both heating the membrane and detecting its temperature. Both resistors are protected by an inert coating in order to prevent chemical reactions with the gas molecules. Above and underneath the membrane, two cavities are etched into the silicon. The measuring gas diffuses into the cavities and depending on the thermal conductivity of the gas more or less energy is dissipated from the membrane (source) which is held at a higher temperature than the remaining chip (hollow). The quantity measured is the voltage that is needed to keep the temperature of the membrane at the higher level. The micro-mechanical chip is fitted into the housing where the measuring gas flows through. In order to allow a temperature controlled operation mode of the sensor, the temperature of housing of the chip is controlled using two heating elements and a temperature sensor that are attached to the housing. The default value is 60°C (140°F).

4 Installation

4.1 General instructions

The **BCP-H2** is protected by packaging on its way to its application location. This secures it against the usual transport strains. However, before installation, check whether the device has been damaged due to improper transport or improper storage. If the device is damaged in any way, operation without hazards is not possible and the device may not be installed and taken into operation.

Check whether the enclosed materials such as seals and screw-caps are suitable for your process conditions (pressure, temperature, etc.).

The installation should only be performed under supervision by a specialist and in compliance with all applicable work safety rules.

4.2 Mechanical connection

Caution!

The sterile filter is not intended for repelling fluids. Never install the sensor such that fluid can run into the measuring adapter.

Ensure that the sensor head and the measuring adapter are correctly positioned.

Protect the measuring adapter from penetration by liquids.

After installation, check that the pipe connection is gas-tight.

4.2.1 Installation on pipes

The connection to a pipe is made with a 1 ¼" nozzle with an external thread:

1. Place the sealing ring (O-ring 30 x 4 mm, viton, item no. Z-OR-00003) on the nozzle (fig. 2).
2. Place on the sensor (fig. 3).
3. Connect the nozzle and the sensor with the screw cap so that the connection is gas-tight (fig. 4).

Note!

Only use the supplied screw caps. Do not use metal screw caps as they result in thermal contact between the measuring adapter and the pipe and thus violate the technical specifications.

Fig. 2

Fig. 3

Fig. 4

4.2.2 Installation on a Tri-Clamp SMS38 connection

:

1. Place the sealing ring (item no. Z-OR-00013) on the nozzle (fig.5).
2. Place on the sensor (fig. 6).
3. Fix the sensor with the Tri-Clamp on the nozzle (fig. 7).

Fig. 5

Fig. 6

Fig. 7

4.2.3 Installation on a POM flow adapter

To install the sensor on a POM flow adapter:

1. Place the sealing ring (item no. Z-OR-00004) on the nozzle of the flow adapter (fig. 8).
2. Place on the sensor (fig. 9).
3. Connect the flow adapter and the sensor with the screw cap so that the connection is gas-tight (fig. 10).

Fig. 8

Fig. 9

Fig. 10

4.2.4 Installation on a stainless steel flow adapter

To install the sensor on a stainless steel flow adapter:

1. Place the sealing ring (item no. Z-OR-00004) on the stainless steel connection piece (fig. 11).
2. Place on the sensor and put the screws in place (fig. 12).
3. Fasten the 4 screws (item no. Z-XX-00007) so that the connection is gas-tight (fig. 13).

Fig. 11

Fig. 12

Fig. 13

4.3 Electrical connection

4.3.1 General information

Caution!

Read the installation instructions carefully to prevent damage to the device.

Proceed step-by-step.

Only use the original plugs, cables and power adapters.

Never connect or disconnect plugs when the device is connected to the power supply.

The device does not have an on/off switch; it starts operation as soon as it is connected to the power supply.

Improper operation can result in damage to the device.

Note!

New wide temperature sensors in aluminum housing (WT sensors) need a power supply with a voltage of 24V. They will not work with 12V. Please check our individual data sheet.

4.3.2 Version 4 – 20 mA in PA6 housing

To connect the measuring device to the connection cable of the sensor head in the PA6 housing (fig. 14), use the supplied socket and strain relief (fig. 15).

Note!

The numbering of the pins and their assignment refer to the socket when seen from behind (fig. 16).

Remove the insulation of the cables a little as possible to avoid short circuits in the plug housing.

Fig. 14

Fig. 15

PIN 1	V+ = 12 V
PIN 2	GND
PIN 3	RS232_TXD
PIN 4	RS232_RXD
PIN 5	1-point calibration
PIN 6	4–20mA, RL < 250 Ohm
PIN 7	For internal use only
PIN 8	GND



Fig. 16: Plug assignment

1. Connect the 12 V DC power supply to pin 1 of the socket.
2. Connect GND (ground) to pin 2.
3. Connect the measuring device to pin 6 ($R_L < 250$ Ohm) and pin 8 GND (ground).
4. Plug the sensor cable into the socket.

After around 1 hour of heating-up time, the sensor still requires adjusting. During the heating-up time, the sensor displays 2.3 mA. To make the adjustment, expose the sensor for approx. 30 minutes (depending on specification – see datasheet) to CO₂ or N₂.

5. Afterwards, connect pin 5 to pin 8 (GND) for 5 seconds.
6. Screw on the strain relief. The sensor has been adjusted.

4.3.3 RS232 serial version in PA6 housing

1. Connect the sensor to the power supply with the cable supplied.
2. Connect the sensor to a computer using the serial cable.


After around 1 hour of heating-up time, the sensor still requires adjusting. To do this, expose the sensor for approx. 30 minutes (depending on specification – see datasheet) to CO₂ or N₂.

The adjustment itself is performed with the **BACVisSingle** software (see **BACVisSingle** operating manual).

Start the **BACVisSingle** software. You will find further relevant information in the corresponding operating manual.

Fig. 17

4.3.4 Version 4 – 20 mA in aluminum housing

To connect the measuring device to the connection cable of the sensor head in the aluminum housing (fig. 19), use the supplied socket and strain relief (fig. 20).

Note!

The numbering of the pins and their assignment refer to the socket when seen from behind (fig. 21).

Remove the insulation of the cables a little as possible to avoid short circuits in the plug housing.

Fig. 20

Fig. 19

PIN 1	V+ see voltage on your individual data sheet
PIN 2	GND
PIN 3	RS232_TXD
PIN 4	RS232_RXD
PIN 5	1-point calibration
PIN 6	4–20 mA, $R_L < 250 \text{ Ohm}$
PIN 7	For internal use only
PIN 8	GND



Fig. 21: Plug assignment

1. Connect the DC power supply to pin 1 of the socket. See voltage on your individual data sheet.

2. Connect GND (ground) to pin 2.

3. Connect the measuring device to pin 6 ($R_L < 250 \text{ Ohm}$) and pin 8 GND (ground).

4. Plug the sensor cable into the socket.

After around 1 hour of heating-up time, the sensor still requires adjusting. During the heating-up time, the sensor displays 2.3 mA. To make the adjustment, expose the sensor for approx. 30 minutes (depending on specification – see datasheet) to CO2 or N2.

5. Afterwards, connect pin 5 to pin 8 (GND) for 5 seconds.

- Screw on the strain relief. The sensor has been adjusted.

4.3.5 RS232 serial version in aluminum housing

- Connect the sensor to the power supply with the cable supplied.
- Connect the sensor to a computer using the serial cable.

After around 1 hour of heating-up time, the sensor still requires adjusting. To do this, expose the sensor for approx. 30 minutes (depending on specification – see datasheet) to CO₂ or N₂.

The adjustment itself is performed with the **BACVisSingle** software (see **BACVisSingle** operating manual).

Start the **BACVisSingle** software. You will find all of the further relevant information in the corresponding operating manual.

Fig. 22

1 = +12V or + 24V, please see you data sheet
2 = 0 V
3 = RS232_RxD
4 = RS232_TxD
5 = RS232_GND = PE
PE = ground

Fig. 23: Plug assignment

4.3.6 Connection via BACCom12

The **BACCom12** connection box is an electronic multiplexer with an integrated pressure sensor. It facilitates the connection of up to 12 sensor heads.

Communication with a PC can be switched between RS232 or Ethernet.

The individual connections are explained in the following table:

	Designation	Description
A	RJ45	RJ45 socket for connecting the sensors
B	LED	Operating display when a voltage is present
C	Sub D 9 pin	Data transmission to the PC
D	Switch	Switches between RS232 and Ethernet
F	RJ45	Ethernet connection
G	Power socket	12 V 3.75 A, only use the supplied power adapter
H	Box reset	Resets the box; does not effect the sensors
K	M8 4 pin socket	4-pin connection sockets A–D for additional boxes

Fig. 24: Front of the **BACCom12**

Fig. 25: Connections on the **BACCom12**

Caution!

To prevent damage to the device, only use the supplied power adapter and the supplied cable.

Never disconnect or connect the connection plugs on the sensor heads when the **BACCom12** is switched on.

1. Connect all sensor heads with the **BACCom12**.
2. Connect the supplied power adapter to the power socket **G**.
3. Plug the power plug of the power adapter into the socket.

After a heating-up time of approx. 1 hour, the measuring system is ready for operation.

4. Connect the **BACCom12** to the PC or network via the Ethernet port **E**, or connect the **BACCom12** via the RS232 output **C** with the **supplied** cable to the serial interface of the computer.
5. Select the corresponding interface with the switch **D**.

After around 1 hour of heating-up time, the sensor still requires adjusting. To do this, expose the sensor for approx. 30 minutes (depending on specification – see datasheet) to CO₂ or N₂.

Adjustment of the sensors is performed with the **BACVis** software. Start the corresponding software. You will find all further information in the software instructions.

After initial commissioning, the measuring device can remain switched on constantly, meaning that the heating-up time is not required before every measuring.

4.4 Cross-sensitivity to other gases

Due to the thermal conductivity measuring principle of the BCP-H2 each sensor has a specific cross-sensitivity to other gases. The value of the deviation depends on the specifications of the used sensor and the gas mixture applied. Every H2 sensor has been pre-calibrated at the factory for a specific dry gas mixture. Exact specifications before the first factory calibration are therefore decisive for the subsequent measurement results. Due to this a BCP-H2 cannot be switched to a different process with other gas components. Please contact us if you want to use it in different gas mixtures.

For detailed information on the cross-sensitivity and the effect and your individual sensor, please ask us directly.

5 Maintenance

We recommend sending the device to BlueSens for annual maintenance, checking and calibration of the sensors.

5.1 1-point calibration

Once monthly, or after each connection and disconnection of the sensor head and measuring adapter, the sensor head must be exposed for approx. 30 minutes (depending on specifications, see datasheet) to CO₂ or N₂.

Afterwards, connect pin 5 to pin 8 on the connection cable for 5 seconds, or, if present on the sensor, press the blue button for 5 seconds (fig. 26).

For the serial version, the adjustment can be made using the **BACVis** software.

Fig. 26

5.2 Recalibration

The sensor should be sent back to the manufacturer or an authorized dealer for annual recalibration.

You can get further information for our annual maintenance service Blue4Care incl. extension of the warranty up to 6 years on our homepage:

<http://www.bluesens.de>

and then in the area:

Service→Blue4Care

5.3 Filter change – coarse filter

5.3.1 Removing the filter cover

Hold the sensor head and screw off the lower cover counterclockwise using the clamping ring tongs (fig. 27). If your filter-cover is already equipped with notches (fig. 28) you can use a big coin (like a 2-Euro-coin or a quarter) to remove the cover, too (fig. 29).

Fig. 27

5.3.2 Changing the filter and seals (Z-XX-00052)

1. Remove the filter from the recess (fig. 30).
2. Insert a new filter.
3. Check the seals for damage and replace as required.
4. Screw on the lower cover using the clamping ring tongs respectively using the coin.

Fig. 28


Fig. 29


Fig. 30

6 Appendix

6.1 Calibration table

Complete calibration can only be conducted by BlueSens. Monthly 1-point calibration can be performed as described in chapter 5.1. The adjustment must also be made each time the sensor head and the measuring adapter are disconnected from each other.

Fill out the table below when this is performed.

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We hereby declare that our product, type:

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complies with the following relevant provisions:

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Low voltage guideline (72/23/EEC and 93/68/EEC) is not applicable as no voltage higher than 24 V is used.

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Applied harmonized standards, in particular:

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