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Operating Manual for the **BlueVCount** gas volume measurement device for measuring volume flows

1.	About this document.....	5
1.1.	Purpose.....	5
1.2.	Target group.....	5
1.3.	Symbols and signal words.....	6
2.	For your safety.....	7
2.1.	In general.....	7
2.2.	Authorised personnel.....	7
2.3.	Proper and intended usage.....	7
2.4.	Storage and transport.....	8
2.5.	Warning against misuse.....	8
2.6.	General safety instructions	9
2.7.	CE and FCC compliance.....	10
3.	Product description.....	11
3.1.	Structure.....	12
3.2.	Measuring volume.....	12
3.2.1.	Measuring principle.....	12
3.2.2.	Converting to a standardized gas.....	12
3.3.	Accessories.....	14
3.3.1.	Included with the delivery of the starter set (SET-STARTBLUECOUNT).....	14
3.3.2.	Additionally required accessories (SET-BLUEVCOUNT).....	15
3.3.3.	Optional accessories	16
4.	Installation.....	17
4.1.	General information.....	17
4.2.	Unpacking and installation	17
4.2.1.	Material.....	19
4.2.2.	Assembly	20
4.2.3.	H ₂ S absorber.....	21
4.2.4.	CO ₂ absorber NaOH	22
4.2.5.	Installation on fermentation vessels.....	23
4.3.	Electrical connection	24
4.3.1.	General information.....	24
4.3.2.	Connection via Modbus RTU.....	25
5.	Operations.....	26
5.1.	At the measurement device	26
5.2.	Using Modbus	27
5.2.1.	Read out measured values.....	27
5.2.2.	Setting values	27
6.	Troubleshooting	28

7.	Maintenance.....	30
8.	Appendix.....	31
8.1.	Technical drawing.....	31
8.2.	Technical data.....	32
8.3.	Warranty	33
8.4.	Service and support.....	33
8.5.	Recycling and disposal.....	33
8.6.	Legal notice.....	34
9.	CE and FCC certificates.....	35
9.1.	CE CERTIFICATE.....	35
9.2.	FCC CERTIFICATE.....	36

1. About this document

1.1. Purpose

This operating manual provides you with all the information you need to quickly commission and safely operate the **BlueVCount**. Please read this operating manual before you start the initial commissioning. Keep this operating manual in a safe place for future reference.

1.2. Target group

This operating manual is intended for trained laboratory personnel (BTA, CTA, etc.). The contents of this manual must be made accessible to the technical staff and implemented by them.

The technical staff must be instructed in and follow the GLP "Good Laboratory Practice" guidelines.

The skilled personnel must have received a safety briefing on the applicable safety regulations.

1.3. Symbols and signal words

This symbol indicates a situation which could result in danger. It is used in combination with the following signal words:

Danger!

A hazard with a high degree of risk.

Failure to avoid this risk will result in death or serious injury.

Warning!

A hazard with a medium degree of risk.

Failure to avoid this risk can result in death or serious injury.

Caution!

A hazard with a low degree of risk.

Failure to comply may result in a minor or moderate injury.

Notice!

Non-observance of notices may result in material damage and impair the function of the product.

This symbol indicates a danger of corrosive substances.

This symbol indicates a danger of poisonous substances.

This symbol indicates a requirement to take an action.

This symbol indicates that the operating manual must be observed.

Notice

This symbol indicates helpful, additional information.

Reference

This symbol indicates additional information referenced outside this document.

List

These list points indicate a list in which the order of the items has no particular significance.

1.

Sequence of steps

These numbers indicate steps that are to be carried out in a certain order.

2. For your safety

2.1. In general

The **BlueVCount** was inspected before it left the factory and was in an operationally ready condition.

Please read this operating manual carefully before installing and commissioning the **BlueVCount**. The operating manual contains safety instructions and warnings that must be followed to ensure safe operation.

Please keep this operating manual in a safe place for future use.

The **BlueVCount** must never be operated under conditions that do not conform to its specifications and information specified on the ratings plate.

Maintenance and repair may only be carried out by properly trained, competent personnel who are familiar with the associated risks and warranty provisions.

2.2. Authorised personnel

The operations described in this operating manual may only be carried out by the trained technical staff who have been authorised by the facility operator.

The technical staff must be instructed in and follow the GLP "Good Laboratory Practice" guidelines.

The technical staff must have received safety training for the applicable accident prevention and safety regulations, and for the gases and liquids being used (Hazardous Substances Ordinance).

For safety and warranty reasons, only personnel from **BlueSens gas sensor GmbH** may intervene in the **BlueVCount** beyond the normal handling steps described in this manual.

2.3. Proper and intended usage

The **BlueVCount** is a measurement device used for measuring volume flows in the specified flow range under the conditions described in the technical data. It is used to monitor metabolism in biological processes and non-biological reactions.

The **BlueVCount** can be used in fermentations, degradation/substrate comparison studies, wastewater treatment studies, bioethanol production, etc. It is particularly useful for measuring gas flows in biogas processes or for residual gas analysis.

Danger!

The **BlueVCount** has no ATEX approval.

There is a risk of serious or fatal injuries.

➔ The **BlueVCount** may only be used in well-ventilated rooms.

2.4. Storage and transport

When it is not in use, store the **BlueVCount** between 0 and 60 °C. The relative humidity should not exceed 75% (non-condensing). Note that dust can clog the **BlueVCount**. Store it in a closed cabinet or box.

Protect the **BlueVCount** from knocks and vibration during transportation

2.5. Warning against misuse

- The **BlueVCount** must not be used as a safety component for monitoring gases in plants.
- It may not be used as a gas detector.
- It may not be used in explosive zones.
- Any use other than the intended usage described in this manual is considered improper use. The manufacturer is not liable for any damage resulting from such usage. The risk lies solely with the operator.

2.6. General safety instructions

The **BlueVCount** may pose application-specific hazards if it is used improperly.

Observe the following general safety instructions:

- Only operate the **BlueVCount** when it is in perfect technical condition.
- When it is not in use, store the **BlueVCount** between 0 and 60 °C. The relative humidity should not exceed 75 % (non-condensing).
Note that dust can clog the **BlueVCount**. Store it in a closed cabinet or box.
- Use a damp cloth to clean the surface of the **BlueVCount**. You can also wipe off the **BlueVCount** using a mild disinfectant.
- If liquid enters the housing, clean the **BlueVCount** immediately.
- Protect the **BlueVCount** from knocks and vibration during transportation.

Danger!

Danger due to possible gas leakage.

There is a risk of severe to fatal injuries due to explosions and poisoning.

- ➔ After installation or adjustment, all connections must be checked for leaks.

Notice!

Please pay attention to the further safety and warning notices in all the following chapters.

Notice

Personnel working with the CO₂ absorber (NaOH) must be protected with suitable protective clothing, namely:

Protective gown,

Safety glasses,

Protective gloves

2.7. CE and FCC compliance

The **BlueVCount** is compliant with the EMC Directive (2014/30/EU) and the RoHS Directive (2011/65/EU).

The Low Voltage Directive (2014/35/EU) does not apply because no voltage higher than 24 V is used.

The CE and FCC certificates can be found in the chapter 9.

3. Product description

The **BlueVCount** is used to measure volume flows.

Overview of the **BlueVCount**:


Figure 1: Front view of the **BlueVCount**


Figure 2: Rear view of the **BlueVCount**

3.1. Structure

The **BlueVCount** mainly consists of the volume measurement unit. This is a cylinder fitted with a magnetic piston that is both sealed and lubricated with oil.

As the magnetic cylinder passes the side-mounted sensors, an electrical pulse is generated and counted. By calibrating the cylinder volume, the flow-through volume can be measured.

The **BlueVCount** itself has a display on which the measured values of the ambient conditions for the pressure and temperature as well as the volume currently measured can be read.

RS485 Modbus is used for communications. The **BlueVCount** can be reset to zero using its reset switch. The Modbus In and Out connections enable several devices to be connected directly to one another ("daisy chained"), so that only one USB RS485 cable is required for connecting to a PC.

The gas supply hoses have an inner diameter of 4 mm and an outer diameter of 6 mm.

3.2. Measuring volume

3.2.1. Measuring principle

The **BlueVCount** works on the principle of a displacement counter. The force on the piston in the guide tube caused by a gas flow acts against the weight of the piston. As a result, a pressure of approximately 3 mbar builds up on the side of the gas source. This gas pressure raises the piston until a magnetic sensor activates the valve and opens the gas outlet. The volume of the gas that has passed through can be calculated from the number of executed piston movements.

3.2.2. Converting to a standardized gas

The **BlueVCount** has a micro-controller which is used to evaluate the magnetic sensors and the valve operations. The micro-controller also converts the volume flow to the standard conditions of 1013.25 mbar and 0 °C. An air pressure sensor is used to detect the reference pressure for the current air pressure at the site location. The temperature at the installation site is measured in parallel. At this temperature, the water vapour content in the gas can be determined based on the water vapour saturation of the incoming gas.

The volume flow of dry gas under standard conditions can be calculated based on these values. The pressure of 3 mbar or higher required for the piston stroke is taken into account in the conversion. On the outlet side of the measurement device, the gas can be discharged without pressure or collected in bags. The **BlueVCount** is capable of forwarding pressure built-up in the gas source to the gas outlet. This pressure increase is not detected by the air pressure sensor. Therefore, the back pressure at the output cannot be taken into account when converting to standard conditions. If there are devices downstream from the **BlueVCount** that build up a back pressure, this back pressure must also be considered when determining the quantity. For example, a downstream drum gas meter will tend to display more gas volume than the **BlueVCount**.

Notice!

Danger from use at excessively high or low temperatures.
There is a risk of material damage and functional impairment.

- ➔ Do not use the **BlueVCount** at temperatures outside the specified operating limits.

Notice!

Danger from using the **BlueVCount** above the specified measuring limits.

There is a risk of material damage and functional impairment.

- ➔ Do not use the **BlueVCount** for volume flows outside the specified measuring limits.

Notice!

Danger from the use of gases containing H₂S.

There is a risk of material damage and functional impairment.

- ➔ Always use the H₂S absorber (Z-XX-00114) for gases which contain H₂S.

Danger!

Danger due to possible explosions.

There is a risk of serious or fatal injuries.

- ➔ Some gases to be measured, such as CH₄, are potentially flammable. Since the **BlueVCount** is not suitable for operation in explosive atmospheres, it must not be operated in such atmospheres.

Danger!

Danger from possible explosions or fires.

There is a risk of serious or fatal injuries.

- ➔ Never use the **BlueVCount** in areas where there is a risk of fire or explosion.

3.3. Accessories

The following accessories are required for the initial commissioning of the **BlueVCount**:

3.3.1. Included with the delivery of the starter set (SET-STARTBLUEVCOUNT)

	Article number	Description	Figure
1.	Z-FM-00031	BlueVCount	 Figure 3
2.	Z-XX-00024	Cold trap 100 ml, incl. GL 45 screw cap (see assembly instructions in section 4.2)	 Figure 4

	Article number	Description	Figure
3.	Z-XX-00114	100 g H ₂ S absorber granulate in a bag	 Figure 5
4.	Z-XX-00109	Polyurethane hose (PUR), 6 x 4 mm, black, 3 m	 Figure 6
5.	Z-BA-MULTI-BVCOUNT	BlueVCount operating manual (this document)	Without figure

Notice!

To ensure that the initial commissioning of the **BlueVCount** is successful, additional accessories are required, which are not included in the scope of delivery. This is described in the following **section 3.3.2 Additionally required accessories (SET-BLUEVCOUNT)**.

3.3.2. Additionally required accessories (SET-BLUEVCOUNT)

	Article number	Description	Figure
1.	Z-KA-00055	An adapter cable for Modbus communications and for the power supply	 Figure 7
2.	Z-KA-00048	USB to RS485 serial port adapter with power supply connection	 Figure 8

	Article number	Description	Figure
3.	Z-NT-00011	24 V DC power supply	 Figure 9
4.	Z-KA-00056	Terminating resistor for Modbus	 Figure 10

Notice!

Optional accessories are available for the **BlueVCount**. This is described in the following **section 3.3.3 Optional accessories**.

3.3.3. Optional accessories

The following accessories are optionally available:

	Article number	Description	Figure
1.	Z-KA-00067	Extension cable, 110 cm	 Figure 11
2.	Z-XX-00121	Set of CO ₂ absorber bottles, without NaOH	 Figure 12

	Article number	Description	Figure
3.	Z-XX-00114	100 g H ₂ S absorber granulate in a bag	 Figure 13

4. Installation

4.1. General information

The **BlueVCount** is protected by special packaging when it is transported to the installation site. (Keep this packaging for maintenance purposes.) This protects the **BlueVCount** from normal transport stresses. Nonetheless, before installation, you should check whether the **BlueVCount** has been damaged by improper transport or storage.

If the **BlueVCount** is damaged in any way, safe operations are no longer possible and the **BlueVCount** must not be installed and put into operation.

Check whether the enclosed materials (such as the seals) are suitable for your process conditions (for the on-site pressure, temperature, etc.).

The installation should only be carried out under the supervision of a professional, in compliance with all applicable work safety regulations.

Make sure that the specified operating conditions are maintained at all times.

4.2. Unpacking and installation

All parts have been securely packaged. Check them for damage while unpacking.

Notice!

Danger due to liquids or dirt particles entering the measuring cylinder.

There is a risk of material damage and functional impairment.

- Make sure that no liquids or dirt particles get into the measuring cylinder of the **BlueVCount**.

Danger!

Danger due to gas leakage.

There is a risk of serious or fatal injuries.

- After installation, check that there are no leaks at the hose connections

Notice!

After the installation, check the gas tightness of the entire system. A leak will result in incorrect measurements.

Notice!

Danger of contamination of the measuring cylinder.

There is a risk of material damage and functional impairment.

- ➔ To avoid contamination of the measuring cylinder, a cold trap must be connected between the fermentation vessel and the **BlueVCount**.

Notice!

Risk of damage to or destruction of the BlueVCount if the following instructions are not followed:

- Observe the **BlueVCount** safety instructions.
- Observe the specifications on your data sheet.
- If the **BlueVCount** has been stored in a cold environment, wait several hours for it to reach room temperature.

Notice!

Danger of decomposition of the sealant of the measuring cylinder.

There is a risk of material damage and functional impairment.

- ➔ In the case of gases containing H₂S, the gas must be purified using an H₂S-absorber granulate.

Notice!

An NaOH solution can be used to absorb the CO₂. However, this is not required to safely use the **BlueVCount**.

4.2.1. Material

The **BlueVCount** comes with a 100 ml glass bottle. This serves to collect the condensate (cold trap) or additionally as an absorber bottle for H₂S. For this, the absorber granulate (not shown) must be added to the bottle (see section 4.2.3).

Letter	Description
A	100 ml glass bottle
B	Screw cap GL 45, with 34 mm hole
C	Screw cap GL 14, with 9 mm hole
D	Locking sleeve
E	Connecting piece for GL 14 screw cap
F	O-ring for connecting piece (E)
G	O-ring for GL 45 bottle
H	Adapter for GL 45, to 4x6 mm hoses
I	Hose (PU) 4mm inner diam. and 6 mm outer diam.

Figure 14

4.2.2. Assembly

1. Place the O-ring (G) on the bottle (A), Figure 15.

Figure 15

2. Then place the hose adapter (H) on the O-ring, Figure 16.

Figure 16

3. Screw the GL45 screw cap (B) on tightly, Figure 17

Figure 17

The hose adapters are marked either **I** or **O**,
Figure 18.

I (beige) means:
"In - gas into the bottle"
(from the fermentation vessel)

O (light grey) means.
"Out - Out of the bottle", towards the
BlueVCount.

Figure 18

4. Tighten the knurled nut and the clamping ring of the hose adapter "I" over the hose (I), Figure 19.

Figure 19

5. Attach the hose to the hose adapter "I", Figure 20.

Figure 20

6. Press the clamping ring into the bulk-head fitting, Figure 21.

Figure 21

7. Tighten the knurled nut, Figure 22.

Figure 22

The other end of the hose must now be connected to the fermenter.

The normally black hose has been replaced by a white hose in the following illustrations for better clarity.

8. Place the red screw cap GL 14 (C) and the locking sleeve (D) on the hose, Figure 23.
9. Attach the nozzle (E) (with O-ring (F)) to the hose, Figure 23.

Figure 23

10. Slide the locking sleeve onto the connecting piece, Figure 24.

Figure 24

11. Attach it to the GL 14 connection of your fermenter, Figure 25.

Figure 25

4.2.3. H₂S absorber

If the absorber is used for H₂S, we recommend that you attach a short hose (10 cm) to the lower part of the input adapter, which is located in the bottle with absorber, Figure 26.

When you insert the tube into the filled jar, you must shake the jar gently, Figure 27.

Figure 26

Figure 27

4.2.4. CO₂ absorber NaOH

Danger!

Danger when using NaOH.

There is a risk of severe skin or eye damage due to chemical burns.

- NaOH may only be used in compliance with the applicable safety precautions.

Personnel working with NaOH must be protected by suitable protective clothing, namely with:

Protective gown

Safety glasses

Protective gloves

The NaOH solution can be used as a CO₂ absorber.

The gas bubbles through the solution and the CO₂ is completely removed with 100 ml 3M NaOH.

To prevent backflow into the fermenter, a cold trap is used as an equalization tank.

Figure 28

Figure 29

Figure 30

1. Prepare two cold traps with 10 cm long hoses touching the bottoms of the glass bottles.

Use a sintered PE silencer as a frit in one of the cold traps.

Figure 29/Figure 30

2. Use a 20 cm PU hose to connect the two hose adapters together (**inlet/in** to **inlet/in**), Figure 31.

This prevents the solution from flowing back into the fermenter when there is negative pressure.

Figure 31

4.2.5. Installation on fermentation vessels

1. Connect the fermentation vessel (a) to the cold trap (c) at the GL 14 hose connection (b) according to the instructions in section 4.2.2.
2. Then connect the hose adapter "O" of the cold trap to the **BlueVCount** (d) using the hose.

Make sure that the hoses are pushed as far as possible onto the hose adapters and then tightened firmly with the knurled nuts.

3. Connect the output of the **BlueVCount** to an exhaust line.

Make sure that this does not create any back pressure.


Figure 32: Example installation for a bio-methane potential (BMP) test system

Figure 33: Example installation for a BMP test station with CO₂ absorber and H₂S absorber

Notice!

Check the intermediate cold trap regularly for any liquid and contamination.

The cold trap must be emptied and cleaned (at the latest) after a completed measurement. If H₂S is generated in your process, the cold trap must be filled with H₂S-absorbing granulate to protect the **BlueVCount**.

To remove the CO₂ contained in the gas, another vessel with CO₂-absorbing NaOH solution can be installed in before the cold trap (see 4.2.4). In this case, another 100 ml cold trap can serve as a backflow protection so that no NaOH can flow back into the fermentation vessel if the flow stops.

Notice!

After installation, check that the screw connection is tight and that there are no gas leaks!

4.3. Electrical connection

4.3.1. General information

Notice!

Risk of damage to the **BlueVCount** due to incorrect installation.

There is a risk of material damage and functional impairment.

- Read the installation manual carefully. Proceed step by step.

Notice!

Risk of damage to the **BlueVCount** due to improper operation/handling.

There is a risk of material damage and functional impairment.

- Follow the instructions for operation and use of the **BlueVCount** as described in this operating manual.

Notice!

Using parts that are not approved for the **BlueVCount** may damage the **BlueVCount**.

There is a risk of material damage and functional impairment.

- Only use the original plugs, cables and power supply units (refer to section 3.3).

Notice!

Risk of damage to the **BlueVCount** due to unintended operations.

There is a risk of material damage and functional impairment.

- The **BlueVCount** does not have an on/off switch; it starts operating as soon as it is connected to the power supply. Make sure you have taken all the necessary safety precautions when connecting the **BlueVCount** to the power supply.

Notice!

Risk of damage to the **BlueVCount** due to a short circuit.

There is a risk of material damage and functional impairment.

- Never plug in or unplug the connectors when the **BlueVCount** is connected to the power supply.

4.3.2. Connection via Modbus RTU

The following items are required to connect the **BlueVCount**:

- a **BlueVCount** adapter cable
- a BlueSens USB-to-RS485 cable
- a 24V power supply
- a terminating resistor.

1. Before connecting the power supply, connect the adapter cable to the RS485 cable and the Modbus IN plug of the **BlueVCount**.

The terminating resistor should then be plugged into the Modbus OUT socket on the other side of the **BlueVCount**.

2. If further measuring devices must be connected, connect them with their Modbus IN plug to the Modbus OUT socket of the previous **BlueVCount**.

Then connect the terminating resistor to the last **BlueVCount**.

Figure 34

Figure 35

Figure 36

3. After all the desired devices have been connected, plug in the power supply unit. This supplies all the devices with voltage.

Figure 37

4. Plug the USB adapter into a free USB port on your computer.

Notice!

The USB adapter does not replace the power supply!

5. After switching on the computer, a COM port should be created automatically in Windows. Some versions of Windows require a driver. Check the USB flash drive (if BlueVis was purchased) or visit the FTDI website to download the latest RS485 driver.

6. If more than one **BlueVCount** is connected to the same bus, the corresponding ModBus IDs must all be unique. The Modbus ID of each **BlueVCount** can be changed using the Modbus (refer to the BlueVCount Communication Manual).

5. Operations

5.1. At the measurement device

The **BlueVCount** gas measurement devices have 2 electrical connections and a push button.

When the device is started, the sensor ID, the current Modbus address and the firmware version (here 1.4.30) are shown on the display (Figure 38 Abb. 38). After approx. 1 minute, the text "BlueSens GmbH" is replaced by the current counter reading (Figure 39).

To reset the unit's counters to zero, press the button once briefly. The display shows "reset^(0)" (Figure 40). By releasing the button (Figure 41) and pressing it down again for more than 1 second, the counter value can also be reset (Figure 42). If the button is not pressed again, the total counter value remains at its previous value (Figure 43).

Since several **BlueVCount** can be operated with one cable set and power supply unit, these must be ordered in the number required for your application. Only original parts may be used; otherwise, the device could be damaged.

After 3 ½ minutes, the display switches to energy-saving mode and only the current counter reading "Total" is displayed. To prevent the display text from burning in, the output regularly moves to another line. Briefly tap the button on the back to display all lines of the output.

RS485 Modbus is used for the communications with the software.


Figure 38

Figure 39
Before resetting

Figure 40

Figure 41

Figure 42

Figure 43

5.2. Using Modbus

Available functions:

- Read holding register (function 3)
- Describe a single register (function 6)
- Writing to multiple registers (function 16)

5.2.1. Read out measured values

The following measured values can be read out:

Register name	Ad-dress	Register	Variable type	Description
Normalized volume [ml]	1	2	Float DCBA	3 decimal places
Ambient pressure [bar absolute]	5	2	Float DCBA	3 decimal places
Temperature [°C]	9	2	Float DCBA	1 decimal place
Sensor ID	13	1	16-bit integer	1 – 65535

5.2.2. Setting values

The following values can be described:

Register name	Ad-dress	Register	Variable type	Description
Set Modbus ID	6656	1	16-bit integer	Range: 1 – 243
Delete total	8000	2	32-bit integer	0

A more detailed description of the Modbus registers can be found in the **BlueVCount** Z-BA-CM-BlueVCount communication manual.

To set up a **BlueVCount** in BlueVis, please use the quick guide for connecting the **BlueVCount** gas volume measurement device to BlueVIS.

6. Troubleshooting

Problem	Possible cause	Solution
The measured values deviate from the expected value.	The flow rate is outside the measuring range.	The BlueVCount has a measuring range due to its design and calibration. Volume flows above this range lead to measurement errors. The measuring range for your device can be found on the ratings plate and in the enclosed data sheet.
	The magnetic piston is no longer moving.	Check that the magnetic piston moves when the BlueVCount is turned upside down. If it does not, contact the service department of BlueSens gas sensor GmbH .
	The NaOH solution is used up (only in connection with CO ₂ absorber).	Replace the used NaOH solution with fresh NaOH solution.
No digital communication	No power supply	Connect the power supply unit.

Problem	Possible cause	Solution
No digital communication	Wrong cable	Verify that the cable harness is complete and properly connected (also refer to section 4.3 Electrical connection).
	The wrong parameters were selected.	Check whether a Modbus RTU connection has been established with the following parameters: Baud rate: 38400 Stop bits: 2 Parity: None

Problem	Possible cause	Solution
No digital communication	Incorrect or duplicate Modbus address/ID	Match the Modbus address of the device on the display with your settings. If two devices have the same address (e.g. factory setting: 1), then each device must be connected individually and assigned a unique address. Refer to section 4.3.2 Connection via Modbus RTU.
	Missing drivers	Check whether the current FTDI driver has been installed (https://ftdichip.com/drivers/).
No volume is measured.	Refer to the problem section "Measured values deviate from the expected value"	Refer to the solutions under problem section "Measured values deviate from the expected value"
Nothing shown on the display	No power supply	Refer to the solutions under problem section "Measured values deviate from the expected value"

Problem	Possible cause	Solution
Nothing shown on the display	Display is defective	Contact the BlueSens service department.


7. Maintenance

To ensure the accuracy of the measured values, the **BlueVCount** should be serviced and calibrated annually.

We recommend sending your **BlueVCount** to **BlueSens gas sensor GmbH** for this.

8. Appendix

8.1. Technical drawing

8.2. Technical data

Application-specific data	
Minimum volume flow	0 ml/h
Maximum volume flow	4,800 ml/h
Mechanical safety mechanism up to	12,000 ml/h
Resolution (for standardized volumes)	< 1 ml
Deviation from the measured value	± 2%

Electrical inputs and outputs	
Power supply	24 VDC
Maximum current consumption	65 mA
Pin assignment for Binder – 5-pin plug	<p>1. +24 V 2. GND 3. RS485_A 4. RS485_B 5. Int. usage</p> 

Operating conditions	
Permissible ambient temperature range	15 – 55 °C
Permissible ambient humidity	< 90 %
Permissible gas temperature	15 – 55 °C
Permissible ambient pressure range:	0.8 to 1.3 bar
Required gas pressure	> 3 mbar

Mechanical specifications	
Width	120 mm (with 5-pin connection)
Height	123 mm
Depth	122 mm (without hoses)
Weight	Approx. 0.81 kg
Protection class	IP64

Electrical input of power supply unit	
Voltage	100 – 240 V ~
Frequency	50 – 60 Hz
Maximum current consumption	0.7 A

Communication interface	
Data format	Modbus RTU
Baud rate	38400
Data bits	8
Stop bits	2
Parity	None
Physical interface	RS485

8.3. Warranty

The warranty period is regulated in the General Terms and Conditions.

Notice

The warranty claim may be voided if this operating manual is not observed, or in case of improper handling (e.g. opening the **BlueVCount**), or if a power supply unit other than the original one is used.

8.4. Service and support

Our qualified customer service is available to you as a partner.

If required, please contact your dealer or us directly:

Phone: +49 2366 4995 567 or e-mail to service@bluesens.de

USA: 847 201 3124 or e-mail to service@bluesens.com

8.5. Recycling and disposal

The **BlueVCount** purchased from us, the accessories and the packaging are all recyclable.

- The **BlueVCount** is packaged to protect it from transport damage. Please take any packaging that you no longer need to your local recycling centre.
- Recycle and dispose of the **BlueVCount**, including its accessories and packaging, in accordance with the national legal requirements applicable.
- The  symbol on the **BlueVCount** means that you are not allowed to dispose of the **BlueVCount** with normal waste within the EU. It must be recycled separately.
- For disposal, we recommend returning the **BlueVCount** and its accessories to **BlueSens gas sensor GmbH**. For this, you will need a **return authorisation number**, which you can obtain from our service team service@bluesens.de together with a valid return address. We guarantee the professional dismantling and environmentally friendly recycling of all components! This service is available free of charge for you.

8.6. Legal notice

BlueVCount operating manual - English

Index Z-BA-MULTI-BVCOUNT

Version V2.00 Rev. 230707 002

© BlueSens gas sensor GmbH

The images and drawings in this operating manual may differ from the original; they are for illustrative purposes only. All information, including the technical information, is subject to change without notice. All images and graphics in this instruction manual:

© BlueSens gas sensor GmbH, Herten.

The specifications in this manual have been carefully checked for accuracy. Nevertheless, **BlueSens gas sensor GmbH** assumes no liability for consequences arising from any errors in the description and illustrations. The General Terms and Conditions from **BlueSens gas sensor GmbH** apply

BlueSens gas sensor GmbH

Snirgelskamp 25

D-45699 Herten, Germany

Phone: +049 (0)2366/4995500

Fax: +049 (0)2366/4995599

E-mail: service@bluesens.de

Home page: www.bluesens.de

Information status: July 2023

9. CE and FCC certificates

The **BlueVCount** complies with the EMC Directive (2014/30/EU) using the harmonised standards EN 55011, EN 61326-1.

The Low Voltage Directive (2014/35/EU) does not apply because no voltage higher than 24 V is used.

9.1. CE CERTIFICATE

EU Declaration of Conformity according to

EU Directive 2014/30/EU (EMC Directive)

EU Directive 2011/65/EU (RoHS Directive)

complies with the essential requirements set out in the Council Directive for the adoption of Member States' legislation on electromagnetic compatibility (2014/30/EU), and the RoHS Directive (2011/65/EU).

The assessment of the product for electromagnetic compatibility was carried out based on the following harmonised EC standards:

EMI EN 55011

EMS DE 61326-1

This declaration applies to all devices manufactured according to the tested prototype.

HERTEN, Germany 7 Nov.
2019

Location and date

Signature

Manufacturer/Representative

We, the **BlueSens gas sensor GmbH**

Sniergelskamp 25, 45699 Herten (Germany)

hereby declare under our sole responsibility that the product:

Gas volume measurement device **BlueVCount**

when used together with our accessories:

Power supply,

Cable

9.2. FCC CERTIFICATE

Declaration of Conformity

We, the **BlueSens gas sensor GmbH**

Snirgelskamp 25, 45699 Herten (Germany)

hereby declare under our sole responsibility that the product:

This declaration applies to all devices manufactured in accordance with
the tested sample.

HERTEN, Germany 7 Nov.
2019

Location and date

Signature

Manufacturer/Representative

Gas volume measurement device **BlueVCount**

when used together with our accessories:

Power supply,

Cable

for which this declaration is based, conforms to the following specifications:

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Wir, die **BlueSens gas sensor GmbH**