Operating Manual (EN)

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Operating Manual

for the **BlueVary**
gas analyser with exchangeable
gas sensor cartridges

1.	About this document.....	5
1.1.	Purpose.....	5
1.2.	Target group.....	5
1.3.	Symbols and signal words.....	6
2.	For your safety.....	7
2.1.	In general.....	7
2.2.	Authorised personnel.....	7
2.3.	Proper and intended usage.....	7
2.4.	Storage and transport.....	8
2.5.	Warning against misuse.....	8
2.6.	General safety instructions.....	9
2.7.	CE and FCC compliance.....	9
3.	Product description.....	10
3.1.	Structure.....	11
3.2.	Measuring principle.....	11
3.2.1.	Pressure and humidity measuring principle.....	11
3.2.2.	Measuring principle of the pure pressure sensor.....	12
3.2.3.	Measuring principle for CO ₂ and CH ₄	13
3.2.4.	O ₂ measuring principle based on ZrO ₂	14
3.2.5.	O ₂ measuring principle based on a galvanic cell.....	15
3.2.6.	H ₂ measuring principle based on thermal conductivity.....	16
3.2.7.	Cartridge compatibility.....	17
3.3.	Accessories.....	19
3.3.1.	Main unit.....	19
3.3.2.	Cartridges.....	21
3.3.3.	Flow adapter.....	22
4.	Installation.....	23
4.1.	General information.....	23
4.2.	Unpacking and installation.....	25
4.2.1.	Mounting the sensor on pipes/hoses.....	28
4.3.	Electrical connection.....	29
4.3.1.	General information.....	29
4.3.2.	Switching on.....	30
4.3.3.	Connection to RS232/USB.....	31
4.3.4.	Connection via Modbus RTU.....	31
4.3.5.	Connection using the analog output.....	32
4.3.6.	Connection via Ethernet Modbus TCP/IP.....	32
4.4.	LCD display.....	34

5.	BlueVary menu.....	35
5.1.	General explanations.....	35
5.2.	Start screen and main menu.....	35
5.3.	Sub-menu for setting up user authorisations.....	36
5.4.	1-point calibration.....	38
5.5.	2-point calibration.....	39
5.6.	Serial interface settings.....	40
5.7.	Modbus settings.....	41
5.7.1.	Modbus RTU settings.....	41
5.7.2.	Modbus TCP settings.....	42
5.7.3.	Modbus RTU Master.....	42
5.8.	Sensor information.....	42
5.9.	Analog settings.....	43
5.10.	Advanced settings.....	45
6.	Troubleshooting.....	46
7.	Maintenance.....	52
7.1.	1-point calibration.....	52
7.2.	2-point calibration.....	53
7.3.	Filter change.....	53
7.3.1.	Removing the sensor cartridges from the flow adapter.....	53
8.	Appendix.....	54
8.1.	Calibration table.....	54
8.2.	Pin assignment of the analog cable art. no. Z-KA-00058.....	54
8.3.	Technical drawing.....	55
8.4.	Technical data.....	56
8.4.1.	BlueVary main unit	56
8.4.2.	BlueVary gas sensor cartridge.....	58
8.5.	Warranty.....	58
8.6.	Service and support.....	58
8.7.	Recycling and disposal.....	59
8.8.	Legal notice.....	59
9.	CE and FCC certificates.....	60
9.1.	CE CERTIFICATE.....	60
9.2.	FCC CERTIFICATE.....	61

1. About this document

1.1. Purpose

This operating manual provides you with all the information you need to quickly commission and safely operate the **BlueVary**.

Please read this operating manual before you start the initial commissioning.

Keep this operating manual in a safe place for future reference.

1.2. Target group

This operating manual is intended for trained laboratory personnel (BTA, CTA, etc.). The contents of this manual must be made accessible to the technical staff and implemented by them.

The technical staff must be instructed in and follow the GLP "Good Laboratory Practice" guidelines.

The technical staff must have received a safety briefing on the applicable safety regulations.

1.3. Symbols and signal words

This symbol indicates a situation which could result in danger. It is used in combination with the following signal words:

Danger!

A hazard with a high degree of risk.

Failure to avoid this risk will result in death or serious injury.

Warning!

A hazard with a medium degree of risk.

Failure to avoid this risk can result in death or serious injury.

Caution!

A hazard with a low degree of risk.

Failure to comply may result in a minor or moderate injury.

Notice!

Non-observance of notices may result in material damage and impair the function of the product.

This symbol indicates a danger of poisonous substances.

This symbol indicates a requirement to take an action.

This symbol indicates that the operating manual must be observed.

Notice

This symbol indicates helpful, additional information.

Reference

This symbol indicates additional information referenced outside this document.

List

These list points indicate a list in which the order of the items has no particular significance.

1.

Sequence of steps

These numbers indicate steps that are to be carried out in a certain order.

2. For your safety

2.1. In general

The **BlueVary** was inspected before it left the factory and was in an operationally ready condition.

Before installing and commissioning the **BlueVary**, read this operating manual carefully. The operating manual contains safety instructions and warnings that must be followed to ensure safe operation.

Please keep this operating manual in a safe place for future use.

The **BlueVary** must never be operated under conditions that do not conform to its specifications and information specified on the ratings plate.

Maintenance and repair may only be carried out by properly trained, competent personnel who are familiar with the associated risks and warranty provisions.

2.2. Authorised personnel

The operations described in this operating manual may only be carried out by the trained technical staff who have been authorised by the facility operator.

The technical staff must be instructed in and follow the GLP "Good Laboratory Practice" guidelines.

The technical staff must have received safety training for the applicable accident prevention and safety regulations, and for the gases being used (Hazardous Substances Ordinance).

For safety and warranty reasons, only personnel from **BlueSens gas sensor GmbH** may intervene in the **BlueVary** beyond the normal handling steps described in this manual.

2.3. Proper and intended usage

The **BlueVary** is a gas sensor for measuring gas concentrations for the specified gases in the specified concentration range under the conditions as described in the technical data. It is used to monitor metabolism in biological processes such as fermentation.

Danger!

The **BlueVary** has no ATEX approval.

There is a risk of serious or fatal injuries.

→ The **BlueVary** may only be used in well ventilated rooms

2.4. Storage and transport

When it is not in use, store the **BlueVary** between 0 and 60 °C. The relative humidity should not exceed 75% (non-condensing). Note that dust can clog the **BlueVary**. Store it in a closed cabinet or box.

Protect the **BlueVary** from knocks and vibration during transportation.

2.5. Warning against misuse

- The **BlueVary** must not be used as a safety component for monitoring gases in plants.
- It may not be used as a gas detector.
- It may not be used in explosive zones.
- Any use other than the intended usage described in this manual is considered improper use. The manufacturer is not liable for any damage resulting from such usage. The risk lies solely with the operator.

2.6. General safety instructions

The **BlueVary** may pose application-specific hazards if it is used improperly.

Observe the following general safety instructions:

- Only operate the **BlueVary** when it is in perfect technical condition.
- When it is not in use, store the **BlueVary** between 0 and 60 °C. The relative humidity should not exceed 75% (non-condensing). Note that dust can clog the **BlueVary**. Store it in a closed cabinet or box.
- Use a damp cloth to clean the surface of the **BlueVary**. You can also wipe off the **BlueVary** using a mild disinfectant.
- If liquid enters the housing, clean the **BlueVary** immediately.
- Protect the **BlueVary** from knocks and vibration during transportation.

Danger!

Danger due to possible gas leakage.

There is a risk of severe to fatal injuries from explosions and poisoning.

- ➔ After installation or adjustment, all connections must be checked for leaks.

Notice

Please pay attention to the further safety and warning notices in all the following chapters.

2.7. CE and FCC compliance

The **BlueVary** is compliant with the EMC Directive (2014/30/EU) and the RoHS Directive (2011/65/EU).

The Low Voltage Directive (2014/35/EU) does not apply because no voltage higher than 24 V is used.

The CE and FCC certificates are listed in chapter 9.

3. Product description

The **BlueVary** is used to measure gas concentrations.


Figure 1: Overview of the **BlueVary**

3.1. Structure

The **BlueVary** sensor consists of a main unit with **three slots for sensor cartridges**. Two different gases can be measured simultaneously. The third slot contains a **pressure and/or humidity sensor cartridge** for automatic compensation. Depending on the application, the gas sensor cartridges are selectable and can be easily replaced. **Yearly maintenance is not required.**


Figure 2: Schematic view of the overall structure

Flow adapters are available for all pipe diameters from 4 mm to 1 1/4". The display of the **BlueVary** shows the status, the measured values from the connected cartridges and other information. This is explained on the following pages.

3.2. Measuring principle

The **BlueVary** analyser may include various measuring components, such as:

- Pressure and humidity measurement
- Pure pressure measurement for aggressive gases
- Infrared sensors for measuring CO₂ and CH₄
- O₂ measurements (based on ZrO₂ or galvanic cell technology)
- H₂ measurement based on the thermal conductivity method

3.2.1. Pressure and humidity measuring principle

A piezo-resistive silicon pressure sensor is used as the pressure sensor. The humidity is measured via a capacitive polymer sensor.

Notice!

Danger from use at excessively high or low temperatures.

There is a risk of material damage and functional impairment.

- ➔ Do not use the **BlueVary** at temperatures outside the specified operating limits.

Notice!

Risk of damage to the humidity or pressure sensor if the **BlueVary** contacts volatile organic compounds.

There is a risk of material damage and functional impairment.

- ➔ Never expose the BlueVary to volatile organic compounds.

Notice!

Risk of damage to the **BlueVary**'s temperature and humidity sensor if exposed to highly concentrated acids or bases or high concentrations of corrosive substances such as NH₃ or H₂O₂.

There is a risk of material damage and functional impairment.

- ➔ Do not expose the **BlueVary** to the substances mentioned above.

3.2.2. Measuring principle of the pure pressure sensor

A stainless steel sensor is used as the pressure sensor. It is immune to aggressive gases, such as those found in biogas.

Notice!

Danger from use at excessively high or low temperatures.

There is a risk of material damage and functional impairment.

- ➔ Do not use the **BlueVary** at temperatures outside the specified operating limits.

3.2.3. Measuring principle for CO₂ and CH₄

The IR cartridge consists of an IR light source, a detector and the reflective measuring cell (Figure 3).

The infrared light beam is reflected by the gas-filled measuring adapter and the light attenuated by the analysed gas is measured by the detector. The translucent sapphire disc prevents the sample atmosphere from escaping and contaminating the optical components. The entire sensor is heated so that no moisture can condense.


Figure 3: Measurement setup of the IR sensor using the example for CO₂

The heating time is approx. 1 hour. This is the case after the initial start-up and each time the **BlueVary** is disconnected from the power supply. The **BlueVary** does not issue any readings during the heating time.

3.2.4. O₂ measuring principle based on ZrO₂

The measuring principle for oxygen for the **BlueVary** is based on an oxygen pump cell that works independently of reference gases. When a voltage is applied to the cell, oxygen ions are pumped from the cathode to the anode.

If the cathode is covered with a gas diffusion barrier, a saturation current forms when the voltage is increased. This is the measure of the oxygen concentration in the atmosphere.


Figure 4: Measuring principle of the O₂ sensor

Notice!

Danger from operating the **BlueVary**'s O₂ gas sensor cartridge based on ZrO₂ in processes without or with insufficient O₂.

Risk of damage to the O₂ sensor element and functional impairment.

→ The ZrO₂-based O₂ gas sensor cartridge always requires a minimum of oxygen so that the O₂ sensor element is not damaged.

The minimum for all concentration ranges is 0.1% vol. % O₂. This minimum level of O₂ must be guaranteed during the entire measurement.

Do not use this cartridge of the **BlueVary** in processes that have no oxygen.

Notice

For measurements that do not contain oxygen in the process, we recommend using an O₂ gas sensor cartridge based on a galvanic cell (O₂ec).

Danger!

Danger from possible explosions or fires.

There is a risk of serious or fatal injuries.

→ Never use the **BlueVary** in areas where there is a risk of fire or explosion.

Notice!

Danger due to silicone components settling on the sensor elements of the **BlueVary**.

There is a risk of material damage and functional impairment.

→ Do not use the **BlueVary** in gases that contain polymer or silicone components.

**Notice!**

Danger from using the **BlueVary** in gases that contain halogens.

There is a risk of damage to the ZrO₂ sensor or dampening pressure sensor.

- ➔ Do not use the **BlueVary** in gases containing halogens such as F, Cl, Br, etc., CFCs or gases containing SOx and H₂S.

**Notice!**

Risk of damage to the sensors due to moisture.

There is a risk of material damage and functional impairment.

- ➔ Never expose the sensor to moisture or water when it is switched off.
(Rel. humidity < 75%)

Notice! Danger when using the **BlueVary** above or below the specified measuring limits.

There is a risk of material damage and functional impairment.

- ➔ Do not use the **BlueVary** in gas concentrations outside the specified measuring limits.

3.2.5. O₂ measuring principle based on a galvanic cell

The oxygen sensor in the cartridge is based on a lead-oxygen battery containing a lead anode, a gold oxygen cathode and a weak acid electrolyte. Oxygen molecules enter the electrochemical cell via a non-porous fluorocarbonide membrane, diffuse into the acid electrolyte and are reduced on the gold electrolytes. The current flowing between the electrodes is proportional to the oxygen concentration in the gas mixture to be measured. The terminal voltages at the thermistor (for temperature compensation) and at the resistor are read as a signal, whereby the change in the output voltages represents the change in the oxygen concentration.


Figure 5: Schematic drawing of the O₂ sensor

3.2.6. H₂ measuring principle based on thermal conductivity

The hydrogen sensor is based on a micro-mechanical silicon chip with a thin membrane. Two thin-film resistors are integrated into the membrane and serve both to heat the membrane as well as for measuring the temperature. Both resistors are protected by an internal coating to prevent chemical reactions with the gas molecules. Two cavities are etched into the silicon above and below the membrane. The sample gas diffuses into the cavities and, depending on the thermal conductivity of the gas, more or less energy is dissipated from the membrane (source), which is kept at a higher temperature than the remaining chip (cavity). The measurement variable is the voltage needed to keep the temperature of the membrane at the higher level. The micro-mechanical chip is installed in the housing through which the sample gas flows. To enable temperature-controlled operation of the sensor, the chip's housing temperature is controlled by two heating elements and a temperature sensor attached to the housing. The default value is 60 °C.

Notice

Cross-sensitivity to other gases!

Due to the thermal conductivity measuring principle of the sensor element, each H₂ cartridge has a specific cross-sensitivity to other gases. The value of the deviation depends on the measurement specifications of the cartridge and the process gas used. For binary gas mixtures, the H₂ cartridge can measure optimally. Each H₂ cartridge has been factory calibrated for a specific dry gas mixture. The exact specification of the operating conditions is therefore crucial for the reliability of the subsequent measurement results, even before the initial calibration at the factory. Therefore, an H₂ cartridge cannot simply be used in another process with different gas components.

Please contact us if you would like to use the H₂ cartridge in other gas mixtures!

For detailed information on cross-sensitivities and possible effects on your individual cartridge, please contact us directly.


Figure 6: Schematic drawing of the H₂ sensor

3.2.7. Cartridge compatibility

Notice!

Risk of damage to cartridges that are not suitable for use in corrosive atmospheres (if this use nevertheless takes place).

There is a risk of material damage and functional impairment.

→ Some gases to be measured, such as biogas, are corrosive. However, some cartridges may not be used with corrosive gases. Therefore, these cartridges may not be used together. If cartridges for corrosive gases are connected together with cartridges that are not suitable for this purpose, the **BlueVary** displays a warning-message during start-up. If the user is sure, they can accept the risk that the life of the unsuitable cartridge will be reduced in this case. For this, the user must enter the admin password. The **BlueVary** will then no longer issue any more error messages for this cartridge.


Figure 7: Error message for incompatible cartridges.

**Danger!**

Danger due to possible explosions.

There is a risk of serious or fatal injuries.

→ Some gases to be measured, such as CH₄, are potentially flammable. Since the **BlueVary** is not suitable for operation in explosive atmospheres, it must not be operated in such atmospheres.

In addition, cartridges for flammable gases must never be used together with cartridges whose sensor element becomes hot, such as the ZrO₂ oxygen cartridge.

**Notice**

When such a cartridge combination is connected to the **BlueVary**, a warning message is shown that cannot be overruled.

Disconnect the **BlueVary** from the power supply immediately!

ATTENTION!

Explosion hazard!

Disconnect from power immediately!

Figure 8: Error message in case of explosion hazard

The table shows the possible cartridge combinations of the **BlueVary**.

Gas compatibility							Required cartridges	
	CO ₂	CO ₂ bio	O ₂ ec	O ₂ (ZrO ₂)	CH ₄	H ₂	Humidity / Pressure	Pressure
CO ₂				✓	✓	✓	✓	✓
CO ₂ bio				✓	✗	✓	✗	✓
O ₂ ec	✓	✓			✓	✓	✓	✓
O ₂ (ZrO ₂)	✓	✗			✗	✗	✓	✓
CH ₄	✓	✓	✓	✗		✓	✗	✓
H ₂	✓	✓	✓	✗	✓		✗	✓

All combinations marked with a ✓ are possible.

All combinations marked with a ✗ may not be used.

The grey fields without marking show a double selection of the same gas. This is also not possible.

3.3. Accessories

The **main unit**, exactly 3 **cartridges** (2 gas sensor cartridges + one pressure and/or humidity sensor cartridge) and a **flow adapter** are required for using the **BlueVary** system.

3.3.1. Main unit

3.3.1.1. Included in the delivery

The delivery of the main unit includes:

	Article number	Description	Figure
1.	K-00015	BlueVary main unit	 Figure 9

	Article number	Description	Figure
2.	Z-KA-00054	Connection cable for BlueVary main unit	 Figure 10
3.	Z-KA-00025	BlueVary connection cable	 Figure 11
4.	Z-KA-00013	Adapter cable RS232 for BlueVary , USB to M12, 1.8 m	 Figure 12

	Article number	Description	Figure
5.	Z-NT-00011	Power supply set (24V power supply unit with international socket adapters)	 Figure 13
6.	Z-XX-00112	Interlock: main unit to cartridge	 Figure 14
7.	Z-BA-MULTI-BVARY	BlueVary operating manual (this document)	Without figure

Notice

Optional accessories are available for the main unit, which are not included in the scope of delivery. This is described in the following section **3.3.1.2 Optional accessories**.

3.3.1.2. Optional accessories

	Article number	Description	Figure
1.	Z-KA-00015	Adapter cable M12 to USB (RS485), 1.8 m	 Figure 15
2.	Z-KA-00057	BlueVary analog cable 2 m, 8 pin	 Figure 16

	Article number	Description	Figure
3.	Z-KA-00058	BlueVary analog cable 5 m, 8 pin	 Figure 17
4.	Z-KA-00059	BlueVary Ethernet cable	 Figure 18
5.		BlueVis software	Without figure

3.3.2. Cartridges

3.3.2.1. Included in the delivery

The scope of delivery of the cartridges includes:

	Article number	Description	Figure
1.		The cartridge itself	 Figure 19
2.	Z-XX-00078	Seal and filter set for BlueVary , single set	 Figure 20

Notice

Optional accessories are available for the cartridges, which are not included in the scope of delivery. This is described in the following section **3.3.2.2 Optional accessories**.

3.3.2.2. Optional accessories

	Article number	Description	Figure
1.	Z-XX-00065	Seal and filter set for BlueVary , 12x 3-pack set	 Figure 21

3.3.3. Flow adapter

3.3.3.1. Included in the delivery

The flow adapter is included in the scope of delivery:

	Article number	Description	Figure
1.		The flow adapter itself	 Figure 22

	Article number	Description	Figure
2.	Z-XX-00079	Retaining plate set for BlueVary flow adapter (metal plates, rubber feet, screws, Torx key)	 Figure 23

Notice

Optional accessories are available for the flow adapters, which are not included in the scope of delivery. This is described in the following section **3.3.3.2 Optional accessories**.

3.3.3.2. Optional accessories

Optional screw-in adapters for the corresponding sizes are available for the flow adapters.

Please contact the sales department of **BlueSens gas sensor GmbH**:

Phone +49 (0)2366 / 4995-500 or e-mail sales@bluesens.de

USA: 847 201 3124 or e-mail to sales@bluesens.com

4. Installation

4.1. General information

The **BlueVary** is protected by special packaging when it is transported to the installation site. (Keep this packaging for maintenance purposes.) This protects the **BlueVary** from normal transport stresses. Nonetheless, before installation, you should check whether the **BlueVary** has been damaged by improper transport or storage. If the **BlueVary** is damaged in any way, safe operations are no longer possible and the **BlueVary** must not be installed and put into operation. Check whether the enclosed materials (such as the seals) are suitable for your process conditions (for the on-site pressure, temperature, etc.).

The installation should only be carried out under the supervision of a professional, in compliance with all applicable work safety regulations.

Notice!

Danger of water entering the measuring cells of the **BlueVary**.

There is a risk of material damage and functional impairment.

- ➔ The integrated filter is not intended for separating liquids (see section 4.2.1). Never install the sensor in such a way that liquid can penetrate.

If water enters the flow adapter, it must be blown out with a sufficient amount of dry air.

Danger!

Danger due to gas leakage.

There is a risk of serious or fatal injuries.

- ➔ After installation, check that there are no leaks at the hose connections.

Notice!

Danger due to contamination of the sensor elements with solvents.

There is a risk of material damage and functional impairment.

- ➔ Do not use solvent-based sealants such as Loctite™ or similar. Please use Teflon or hemp tape as sealing material when needed.

Notice!

Risk of damage to or destruction of the BlueVary if the following instructions are not followed:

- Observe the **BlueVary** safety instructions.
- Never start the **BlueVary** if it has been exposed to high levels of humidity.
- Observe the specifications on your data sheet.
- Do not start the **BlueVary** with an excessive concentration of moisture in the measuring chambers or condensation in the flow adapter.
- Never expose the sensor to water or moisture when it is switched off (rel. humidity < 75%).
- If the **BlueVary** has been stored in a cold environment, wait several hours for it to reach room temperature.

4.2. Unpacking and installation

All parts have been securely packaged. Check them for damage while unpacking.

In order to assemble and mount the **BlueVary**, you will require (Figure 24 and Figure 25) the following:

1. The **BlueVary** main unit (A), (art. no. K-00015)
2. A flow adapter (B), (various types available)
3. Two cartridges for gases and one cartridge for pressure/humidity C), (various types available)

Figure 24: Components of the **BlueVary**

4. Interlock: main unit to cartridge (D), (art. no. Z-XX-0012)
5. Filter set and seal for **BlueVary** (E, F, G), (art. no. Z-XX-00078)
6. Retaining plate set for **BlueVary** flow adapter (Art.-No. Z-XX-00079)

Figure 25: **BlueVary** filter and seal

Mount the **BlueVary** as follows:

1. The **BlueVary** can be mounted vertically (A) or simply on a base (B) (Figure 26).
2. Do not attach the rubber feet if you want to mount the **BlueVary** to a wall. Otherwise, carry out the following steps.
3. Attach the rubber feet to the outer holes of the mounting plate (Figure 27) using the larger screws.


Figure 26: Flow adapter with mounting plates for horizontal mounting or placement on a surface

Figure 27: All parts of the mounting plates / rubber feet already mounted

- Secure the mounting plates to the flow adapter (Figure 28) using the remaining screws.

Figure 28: Fastening the mounting plates to the flow adapter

- Place the seals in the recesses (Figure 29).
- Insert the thin filter into the recesses (Figure 30).
- Insert the thicker filters into the recesses (Figure 30).
- If the filters do not fit into the recesses, you may have to press them in slightly.


Figure 29: **BlueVary** filter and seal

Figure 30: Filters and seals used on the **BlueVary**

Next, place the filters and the seals on the flow adapter.

All necessary filters and seals are included in the filter set for the **BlueVary** (art. no.: Z-XX-00078). A seal and two different filters are provided for each of the three recesses of the flow adapter. There is a thin filter (Figure 25: E), a thicker filter (Figure 25: F) and a seal (Figure 25: G). The thin filter is hydrophobic (water repellent) and the thicker filter is designed to retain impurities and particles. For each recess, there is a corresponding seal, a thin filter and a thicker filter in the filter set (Figure 25).

9. The gas sensor cartridges can now be connected to the flow adapter. To do this, place the first cartridge at an angle of 10 to 20° on the right side of the flow adapter and push it slightly in the direction of the flow adapter until the cartridge noticeably lowers into the recess around the seal (Figure 30). Then, press it against the flow adapter and turn it horizontally (Figure 31).
10. Place the second cartridge to the left of the first cartridge and, like the first cartridge, push it into the recess on the flow adapter and turn it (Figure 32).
11. The third cartridge is placed in the same way onto the last free space on the flow adapter (Figure 33). This cartridge must be the first to be removed when changing cartridges!

Figure 31: The first cartridge must be placed on the right side of the flow adapter.

Figure 32: The second cartridge must be placed to the left of the first cartridge at an angle of 10°.

Figure 33: The third cartridge must also be placed at an angle of 10° from left to right. When removing the cartridges, start from left to right.

12. Then, the **BlueVary** main unit can be placed on the cartridges. For this, the cartridges must be positioned straight on the flow adapter. Then, the **BlueVary** main unit is connected onto the plugs of the cartridges.

Figure 34: Placing the **BlueVary** main unit on the cartridges

13. To connect the **BlueVary** main unit and the cartridges securely, the main unit / cartridge lock mechanism is inserted into the top of the **BlueVary** main unit.

Figure 35: Inserting the lock mechanism for the main unit / cartridge

14. When the main unit / cartridge lock is fully inserted, the curved pin at the end of the interlock is engaged in the notch provided in the **BlueVary** main unit.

Figure 36: Engaged interlock for the main unit / cartridge

4.2.1. Mounting the sensor on pipes/hoses

The sensor has a flow adapter (Figure 24: B) through which the gas to be measured flows. It is available in different sizes from 1/8" to 1 1/4".

Various connectors for pipes and hoses are available for sizes 1/8" to 1/4". This allows the sensor to be connected directly to the exhaust gas line of fermenters.

Filters for particle retention are located between the flow adapter and the actual measuring cell. These filters are not impermeable to water.

Figure 37: Gas flows in the **BlueVary** "from top to bottom" so that liquids can be flushed out.

The sensor should always be mounted so that **no** liquid can penetrate the measuring cells. Never install the **BlueVary** at the lowest point of the line. Make sure that the gas can flow in a vertical direction (from top to bottom) through the flow adapter (Figure 37).

4.3. Electrical connection

4.3.1. General information

Notice!

Risk of damage to the **BlueVary** due to incorrect installation.

There is a risk of material damage and functional impairment.

- Read the installation manual carefully. Proceed step by step.

Notice!

Risk of damage to the **BlueVary** due to unintended operations.

There is a risk of material damage and functional impairment.

- The **BlueVary** does not have an on/off switch; it starts operating as soon as it is connected to the power supply. Make sure you have taken all the necessary safety precautions when connecting the **BlueVary** to the power supply.

Notice!

Using parts that are not approved for the **BlueVary** may damage the **BlueVary**.

There is a risk of material damage and functional impairment.

- Only use the original plugs, cables and power supply units (Figure 38).

Notice!

Risk of damage to the **BlueVary** due to improper operation/handling.

There is a risk of material damage and functional impairment.

- Follow the instructions for operation and use of the **BlueVary** as described in this operating manual.

Figure 38: Original cable and electrical accessories

Figure 39: Port A:
Power supply, RS232, ModBus RTU
Port B: 3 x analog outputs, ext. 1-point
input

4.3.2. Switching on

1. Select the appropriate socket adapter B (Figure 38) that is required in your country for connecting the accompanying mains adapter A (Figure 38) (both art. no. Z-NT-00011).
2. Connect the cable E (Figure 38), art. no. Z-KA-00054, to connection A (Figure 39).
3. Connect the cable E, art. no. Z-KA-00054, and cable C, art. no. Z-KA-00025 (use plugs of both cables).
4. Connect the power supply unit A, art. no. Z-NT-00011, to cable C, art. no. Z-KA-00025 (Figure 40).

Figure 40: Connecting the power supply A to the plug of cable C

4.3.3. Connection to RS232/USB

The following components are required for the electrical connection of the RS232 interface to USB (Figure 41):

Figure 41: Blue RS232 sticker on the USB adapter of cable D

1. Connect cable D, art. no. Z-KA-00013, to the corresponding plug of cable C.
2. Plug the USB adapter into a free USB port on your computer. Make sure that the USB adapter does not replace the power supply!
3. After switching on the computer, a COM port should be created automatically in Windows. Some versions of Windows require a driver. Check the USB flash drive (if BlueVis was purchased) or visit the FTDI website to download the latest RS232 driver.
4. **Perform a 1-point calibration for each gas sensor cartridge before use! See section 7.1.**

4.3.4. Connection via Modbus RTU

The following parts are required for the electrical connection: Cable RS485 to USB, Article number: Z-KA-00015 (not included in the scope of delivery).

Figure 42: Green RS485 sticker on the RS485 USB port Modbus RTU

1. Connect the RS485 cable to the corresponding plug of cable C.
2. Plug the USB adapter into a free USB port on your computer. Make sure that the USB adapter does not replace the power supply!
3. After switching on the computer, a COM port should be created automatically in Windows. Some versions of Windows require a driver. Check the USB flash drive (if BlueVis was purchased) or visit the FTDI website to download the latest RS485 driver.
4. **Perform a 1-point calibration for each gas sensor cartridge before use! See section 7.1.**
5. If more than one **BlueVary** (up to 12) are connected to the same bus, the corresponding ModBus IDs must all be unique. The Modbus ID of each **BlueVary** can be easily changed using the menu. (Refer to section 5.) **The RS485 data bus must be terminated with a terminating resistor (part of the Z-KA-00015) at the end of the data bus.**

4.3.5. Connection using the analog output

For the analog connection of the **BlueVary** to older process control systems, the analog cable Z-KA-00058 is required (not included in the scope of delivery), (for the pin assignment, refer to section 8.2).

1. **Perform a 1-point calibration for each gas sensor cartridge before use! See section 7.1.**
2. The analog outputs must be set via the corresponding menu. The values for offset and slope may be selected as needed; they are not set by the connected cartridges! See section 5.9.

4.3.6. Connection via Ethernet Modbus TCP/IP

The **BlueVary** has Modbus TCP/IP connectivity so that it can easily be used in corporate networks.

1. Switch on the **BlueVary** as described in section 4.3.2.
Connect the Ethernet cable (Figure 45) to connection C (Figure 44).
2. Connect the RJ45 plug to your network. The IP settings for the **BlueVary** are described in chapter 5.
3. **Perform a 1-point calibration for each gas sensor cartridge before use! See section 7.1.**

Figure 43:
Analog cable with open
ends

Figure 44:
Connection C: M12 Ether-
net for ModBus TCP/IP
connection D: For future
usage

Figure 45:
Ethernet cable for Mod-
Bus TCP/IP

Notice!

Risk of damage to or destruction of the *BlueVary* if the following instructions are not followed:

- Observe the ***BlueVary*** safety instructions.
- Never start the ***BlueVary*** if it has been exposed to high levels of humidity.
- Observe the specifications on your data sheet.
- Do not start the ***BlueVary*** with an excessive concentration of moisture in the measuring chambers or condensation in the flow adapter.
- Never expose the sensor to water or moisture when it is switched off (rel. humidity < 75%).
- If the ***BlueVary*** has been stored in a cold environment, wait several hours for it to reach room temperature.

4.4. LCD display

After connecting the **BlueVary** to the power supply (refer to section 4.3), the cartridges are queried. Communication is established between the cartridges and the **BlueVary** main unit:

If no error occurs, the start display appears after a few seconds. The background colour changes to yellow and remains yellow during the heating time of approx. 60 minutes. During this time, no concentration measurements are possible. After heating up, the background colour changes to white*. The **BlueVary** is now ready for use!

Notice

After switching on for the first time and after each cartridge change, the **BlueVary** requests a 1-point calibration. This is signalled by a yellow/green flashing background colour. (1-point calibration: refer to section 7.1).

Figure 46: Display when searching for cartridges

Figure 47: Start display during the heat-up time


Figure 48: Start display after heating up (and 1-point calibration)

Depending on the status of the **BlueVary**, the background can have different colours:

Colour	BlueVary status
White	Normal state, BlueVary ready for operation
Yellow	BlueVary in the heating phase (60 min. after switching on)
Blue	Communication via serial interface
Yellow/ green flashing	1-point calibration of at least one cartridge required
Yellow/ white flashing	Remaining life of at least one cartridge less than 1 month (time to order a new cartridge)
White/ red flashing	Service life of at least one cartridge expired or device error (service required)

* The white background colour may also appear light blue depending on the viewing angle.

5. BlueVary menu

5.1. General explanations

The **BlueVary** has an integrated keyboard and LCD display for displaying values and changing settings. To mark the selected line, a "-" [cursor] is displayed in front of the text of the line (see Figure 49).

The functions of the keyboard are:

- ↑: Moves the cursor up / increases values
- ↓: Moves the cursor down / decreases values
- ←: Moves the cursor to the left
- : Moves the cursor to the right
- OK**: Make input confirmation / selection
- esc**: End of input / back to last level

Figure 49:
Available buttons on the
BlueVary main unit

- To select a menu item, move the cursor to the corresponding line using ↑↓ and press **OK**.
- If there are options, make your selection using ↑↓.
- Enter numbers using ↑↓←→.
- To cancel an entry, press **esc**.
- To apply the values, select the menu item
- *Apply settings*.
- If you want to leave a menu, you can press **esc** at any time.

5.2. Start screen and main menu

The start screen is shown in Figure 50. It displays the gas type and concentration of the connected cartridges, followed by the pressure and absolute humidity. The remaining service life of the connected cartridges is also displayed.

CO₂ : 0.042 vol.%
O₂ : 20.955 vol.%
Pres: 1.011 Bar
Hum : 1.095 vol.%
CO₂ : 339 days left
O₂ : 325 days left
Hum : 701 days left
-Enter = Menu

CO₂ : 0.040 vol.%
O₂ : 20.978 vol.%
Pres: 0.982 Bar
Hum : 1.095 vol.%
CO₂ : 344 days left
O₂ : 352 days left
Hum : 729 days left
-Enter = Menu

Glossary of figures:

Bold lines are for input or selection.

The designated button is printed in **bold** in the text: ↑↓←→ **OK** **esc**

Menu items are printed in *italics* in the text.

Figure 50: Start screen, scheme and image

Press **OK** to go to the main menu:


Figure 51: Start screen and the sub-items of the main menu (schematic)

5.3. Sub-menu for setting up user authorisations

In the sub-menu *Set user level*, user-based access rights can be selected according to the following list:

User level	Access rights	Code
Normal user	Read rights	Unnecessary
Calibration user	Normal user + calibration rights (1-point/ 2-point)	0001
Installation user	Normal user + configuration rights	2020
Admin	Unlimited rights	1818

Notice

After a reset or switch-on, the user level always changes to *Normal user*!

To select a user, first position the cursor in front of the user you want to select using the **↑↓** buttons. Press **OK** to select the user. A code must be entered for any user other than the normal user. Highlight the code line (**↑↓ + OK**) and enter the corresponding code (**↑↓←→ + ESC**).


Figure 52: Schematic of the "User level" sub-menu and display during code entry

To accept the entry, go to *Apply settings* (**↑↓**) and press **OK**.

Depending on whether the code was correct or not, one of the following messages will be displayed:

Figure 53: Possible messages after trying to apply settings in the user menu

Confirm the message with **OK**.

5.4. 1-point calibration

Notice

The 1-point calibration can only be performed when the *calibration user* or *admin* is logged in.

In the *1-point calibration* sub-menu, the cartridges to be calibrated can be selected first. Only one or both cartridges can be selected at the same time. The corresponding gas concentrations required for the 1-point calibration are also displayed.

Figure 54


Figure 55: Schematic of the sub-menu "1-point calibration" and display after selection of both cartridges.

The error 'Deviation over the last 10 minutes...' is issued if the signal has not yet stabilised or deviates too far from the set-point value (e.g. due to a poor last 1-point calibration).


Figure 56: Error message for 1-point calibration

The error 'Not enough values ...' (for calibration) is issued if the **BlueVary** has been disconnected from the power supply or reset within the last 10 minutes. No meaningful 1-point calibration can be carried out.


Figure 57: Error message for 1-point calibration

Notice

For more information on 1-point calibration, see section 7.1

5.5. 2-point calibration

Notice

The 2-point calibration can only be performed when the *calibration user* or *admin* is logged in.

The structure of the *2-point* calibration sub-menu is almost identical to the *1-point* calibration. The first difference is that only one cartridge can be calibrated at a time. The second difference is that the concentration for the 2-point calibration can and must be entered by the user.


Figure 58: Schematic of the "2-point calibration" sub-menu

Notice

If the last 1-point calibration was more than 10 hours ago, a 2-point calibration is not possible. The following error message appears:

1-point too long ago
for 2-point
-Cancel

Figure 59: Error message for 2-point calibration

Notice

The 2-point calibration should only be carried out by experienced users after a perfect 1-point calibration! Please contact the service department of **BlueSens gas sensor GmbH** for more information!

5.6. Serial interface settings


Notice

The serial interface settings can only be changed when the *installation user* or *admin* is logged in.

Notice

The settings displayed correspond to the factory settings as delivered.

The value can only be changed in steps of 10 and corresponds to the seconds of the interval. To stop the periodic output, specify 000 here.

Figure 60: Schematic of the "Serial settings" sub-menu

In the sub-menu *Serial settings*, the settings of the serial interface (RS232) can be changed. The baud rate can be set to 9600, 19200 or 38400 baud. In addition, a periodic output mode can be selected. To activate the periodic output mode, in the line *Send data every:* specify a value greater than 000.

5.7. Modbus settings

Notice

The Modbus settings can only be changed when the installation user or *admin* is logged in.

Notice

The settings displayed correspond to the factory settings in the delivered state.


Figure 61: Schematic of the "Modbus settings" sub-menu

In the *Modbus settings* sub-menu, the settings for the Modbus interface (RTU/RS485 or TCP/IP) can be changed. There are three additional sub-menus (*Modbus RTU settings*, *Modbus TCP settings* and *Modbus RTU master*). There is also the option to select the *Modbus Mode*, set the *Modbus id* and select the *Extended Mode* to be able to communicate with certain Infors fermenter controllers.

Notice

Only set *Extended Mode* if the **BlueVary** is directly connected to an Infors fermenter controller. If in doubt, ask your Infors service representative whether your fermenter requires this setting.

Depending on the selected Modbus mode, the corresponding settings can be changed in the sub-menus (RTU or TCP).

5.7.1. Modbus RTU settings

Notice

The Modbus RTU settings can only be changed when the *installation user* or *admin* is logged in.

Notice

The settings displayed correspond to the factory settings as delivered.


Figure 62: Schematic of the "Modbus RTU settings" sub-menu

In the *Modbus RTU settings* sub-menu, the settings of the Modbus RTU interface (RS485) can be changed. The baud rate can be set to 9600, 19200 or 38400 baud. The parity can be set to "none" (with 2 stop bits) or "even" (with 1 stop bit). In addition, a special Modbus extended mode can be enabled. This allows the **BlueVary** to be connected directly to some Infors fermenters. Please contact your local Infors dealer for further assistance.

5.7.2. Modbus TCP settings

Notice

The Modbus TCP settings can only be changed when the *installation user* or *admin* is logged in.

Notice

The settings displayed correspond to the factory settings as delivered.


Figure 63: Schematic of the "Modbus TCP settings" sub-menu

In the sub-menu *Modbus TCP settings*, the settings for the Modbus TCP interface (TCP/IP) can be changed. The *IP* (address), *Mask* (subnet mask) and *Gate* (gateway) can be set. The *Ph* (MAC address) is only displayed and cannot be changed.

Notice

Please consult your network administrator before changing settings or switching to Modbus TCP mode!

5.7.3. Modbus RTU Master

Notice

For future use!

5.8. Sensor information


Figure 64: Schematic of the "Sensor Info" sub-menu

The sub-menu *Sensor Info* shows additional information about the connected cartridges, the *Main unit (BlueVary main unit)* and the installed *firmware* (version). You can get additional information by selecting the desired cartridge or main unit (the *third channel* is only intended for future usage).

5.9. Analog settings

Notice

The analog settings can only be changed when the *installation user* or *admin* is logged in.

Notice

The settings displayed correspond to the factory settings as delivered.

Notice

During initial installation, use Set Output to 4 mA and Set Output to 20 mA to adapt your PLC or process controller to the **BlueVary**.


Figure 65: Schematic of the "Analog settings" sub-menu

In the *Analog settings* sub-menu, the three analog outputs of the **BlueVary** can be configured. Since the values do not change automatically when cartridges are changed, they may need to be adjusted manually according to the user's needs. The outputs can also be set to 4 mA or 20 mA to control external devices such as a PLC or process controller.

The first two channels correspond to the built-in gas sensor cartridges (here CO₂ and O₂). The input for the *third channel* can be selected (here *humidity*).


Figure 66: Schematic of the sub-menus for the first two channels (top) and the third channel (bottom)

The user can choose to enter the *min.* and *max.* values or the *offset* and *slope* (mA/%). When *min.* or *max.* values are entered, the *offset* and *slope* values are automatically calculated and displayed. Upon entering *offset* or *slope* values, the *min.* and *max.* values are automatically calculated and displayed. For the *third channel*, the *input* can be selected as *humidity* (absolute), *pressure* (absolute) and *external* (for future usage only).

5.10. Advanced settings

Notice

The selection of an external sensor is intended only for future use.

Notice

The external 1-point input (refer to section 7.1) can only be activated if the *installation user* or *admin* is logged in.

Notice

The reset of the 2-point calibrations is intended only for future use.

Notice

Overwriting the memory with a new firmware can only be done when the *admin* is logged in.


Figure 67: Schematic of the "Advanced settings" sub-menu

The **Advanced settings** sub-menu enables the user to access to some additional functions of the **BlueVary**. Under the menu item **Ext. 1-pt** [External 1-point], the connection for external triggering of the 1-point calibration can be activated (see section 7.1). This makes it possible to trigger the 1-point calibration of all gas cartridges of the **BlueVary** simultaneously using the relay of a control unit.

The **Read new firmware** option loads new firmware which enables the **BlueVary** to benefit from new features that will be available in future updates. In order to use this function, the new firmware must be downloaded together with the latest bootloader software from the service page of the **BlueSens gas sensor GmbH** website.

Notice

The menu settings of the **BlueVary** are retained after a firmware update.

Notice

Please contact the **BlueSens gas sensor GmbH** service department for more information!

Notice!

Risk of damage to the **BlueVary** due to handling errors such as disconnecting the power supply during the firmware update

There is a risk of serious damage to the BlueVary or the BlueVary becoming unusable.

→ Do not disconnect the **BlueVary** from the power supply during the firmware update.

6. Troubleshooting

Problem	Possible cause	Solution
Display lights up/flashes red. No measured values or not all values are displayed	No communication with the cartridges is possible	Disconnect the main unit from the power supply and disconnect the cartridges from the main unit. Check that the contacts (in the round socket) are intact and straight. Reassemble the components and test again.
	Insufficient power supply	Check that the power supply unit outputs 24V DC 1A. Always use original accessories from BlueSens gas sensor GmbH (art. no. Z-NT-00011).
	Firmware is not up to date	Newer cartridges require a firmware update. You can get the firmware update in the service page of the BlueSens gas sensor GmbH website (www.bluesens.com).

Problem	Possible cause	Solution
Display lights up/flashes red. No measured values or not all values are displayed	Incompatible cartridges have been combined	Some cartridges must not be combined for safety reasons. You will find a list of the permitted cartridges in this operating manual in the section 3.2.7 cartridge compatibility.
Display flashes alternately green/yellow	At least one cartridge needs a 1-point calibration.	Carry out a 1-point calibration according to the operating manual (see section 7.1 1-point calibration).
Display is orange/yellow and some measured values show NaN.	Sensor heats up	After the initial commissioning, the cartridges need about 60 minutes to reach a stable operating temperature. No measured values are output during this time.
	Insufficient power supply	Check that the power supply unit outputs 24V DC 1A. Always use original accessories from BlueSens gas sensor GmbH (art. No. Z-NT-00011).

Problem	Possible cause	Solution
Display is orange/yellow and some measured values show NaN.	Power supply is unstable	Install a sufficiently capable UPS (uninterruptible power supply) between the socket and the power supply unit.
Measured values deviate from the expected value	Expected value is outside the measuring range.	Cartridges from BlueSens gas sensor GmbH have a predefined measuring range, which is specified on the ratings plate and in the enclosed data sheet. Measured values outside this range are displayed, but are outside the specified accuracy. If your process is regularly out of range, you can simply replace the affected cartridge with a new cartridge that has a different measuring range. Simply contact the sales department of BlueSens gas sensor GmbH .

Problem	Possible cause	Solution
Measured values deviate from the expected value	Missing or incorrect 1-point calibration	Transport stress or a 1-point calibration that has not been carried out properly may cause a strong drift of the signal. Please carry out a new 1-point calibration (refer to section 7.1 1-point calibration).
	Incorrect 2-point calibration	An poorly chosen second point or an improperly performed 2-point calibration can also result in a strong drift. Delete the 2-point calibration (see section 5.5 2-point calibration) and check whether the drift still exists. If it does, please perform a new 1-point calibration.

Problem	Possible cause	Solution
No digital communication	Wrong cable	Please verify the completeness and proper connection of the wiring harness (see section 4.3 Electrical connection) and use the cable with the proper label for your protocol:
		BlueSens serial protocol (RS232): You will need the USB cable with the RS232 label (label on M12 and USB plug, art. no. Z-KA-00013).
		Modbus RTU: Use the USB cable labelled RS485 (label on M12 and USB connector, art. no. Z-KA-00015). You must connect the terminating resistor to the last free socket of the daisy chain.

Problem	Possible cause	Solution
No digital communication	Wrong cable	Modbus TCP: Use the original BlueVary RJ45 connection cable (art. no. Z-KA-00059).
		Wrong parameters selected Compare your settings with those in the BlueVary : Section 5.6. Serial interface settings for RS232, Section 5.7.1 Modbus RTU settings for Modbus RTU and Section 5.7.2 Modbus TCP Settings for Modbus TCP. The IP address of the BlueVary must be in the same network as the receiver and must not be assigned twice. If in doubt, consult your IT department for the appropriate settings.

Problem	Possible cause	Solution
No digital communication	Modbus has wrong operating mode	The BlueVary cannot be addressed simultaneously via Modbus RTU and TCP. Please select your correct operating mode in the Modbus settings.
	Missing drivers	For RS232 or RS485, check whether the latest FTDI driver has been installed (https://ftdichip.com/drivers/).
Improper or no 4-20 mA output	Wrong cable	Please use only original accessories from BlueSens gas sensor GmbH for the connection (article no. Z-KA-00058).
	Wrong pin assignment	Please compare the pin assignment with the table under section 8.28.2 Pin assignment of the analog cable, art. no. Z-KA-00058.

Problem	Possible cause	Solution
Improper or no 4-20 mA output	Incorrect settings	The BlueVary permits its three analog outputs to be freely configured. Therefore, check whether the settings on your control system/software match the values for the min/max, offset and slope. Note that for the third channel, the value that may have been transferred must be selected (see section 5.9 Analog settings).
Missing calibration with the receiver (analog calibration)		During the initial setup, the analog output of the BlueVary must be aligned with the input of your controller. For this, the outputs of the BlueVary can be set to 4 or 20 mA (refer to section 5.9 Analog settings).

Problem	Possible cause	Solution
Error 1-Pt	Outdated firmware	This error message may be caused by an error in the firmware. Please download the current firmware version from the service page of the BlueSens gas sensor GmbH website (www.bluesens.com).
	Impurities/leaks in the supply line for the test gas	The error is triggered by fluctuations in the sensor's signal. This can be caused by leaks or contamination. Especially when working with compressed air, it is essential to check whether the quality meets the necessary requirements.
	Wrong gas	This error can also be triggered by an excessive difference between the measured and the target concentrations. Either the last 1-point was incorrect, or the current gas is incorrect.

Problem	Possible cause	Solution
Wrong gas	Test gas cylinder is empty	Check the pre-pressure of your test gas cylinder and replace it if necessary
The BlueVary does not respond to inputs.	The BlueVary is in an improper operating state.	If the BlueVary does not respond to any input, disconnect the BlueVary from the power source for at least 30 seconds. After that, it should work normally again. If such failures occur frequently, please contact the service department of BlueSens gas sensor GmbH .

Problem	Possible cause	Solution
The BlueVary does not respond to inputs.	Insufficient user level	If the BlueVary responds to inputs but you cannot perform certain actions, your assigned user level is not sufficient. Check in the operating manual to find out which level is necessary for the desired action. Change to this level as described in the manual. You will also find the necessary passwords there.
No display, no illumination after firmware update	Incomplete update, error during update	If the update is incomplete or an error has occurred during the process, the BlueVary remains in update mode. Please reconnect the update software and carry out an update.

7. Maintenance

The **BlueVary** is maintenance-free because of its replaceable gas sensor cartridges. The individual cartridges for the gases have an operating lifespan of 365 days. The remaining lifespan is continually shown on the display of the **BlueVary** main unit. Cartridges for humidity/pressure or pressure have a lifespan of 730 days until they need to be changed. The **BlueVary** main unit is maintenance-free.

In addition to the regular 1-point calibration as described in section 7.1, we recommend replacing the filters depending on the contamination.

7.1. 1-point calibration

The sensor needs a 1-point calibration to maintain measurement accuracy once per month. For this, expose the sensor to the gas specified in the respective sensor cartridge data sheet for at least 30 minutes. The flow rate must be at least 200 ml per minute. If the **BlueVary** has been exposed to high concentrations of CO₂, flushing for several hours is necessary. To conserve calibration gas in this particular case, normal fresh air can be used first and the specified gas only in the last 30 minutes of the purging process. Note that NOT all **BlueVary** sensor cartridges are suitable for 1-point calibration with fresh air. Changes in the gas mixture during 1-point calibration will lead to errors in the measurement. The standard gas concentrations for the **BlueVary** cartridges are listed in the table below. Please also refer to the specific data sheet of your **BlueVary** cartridge. Information on the data sheet has priority over all other information!

Standard concentrations for 1-point calibration:

(To find out which gas your sensor is calibrated for, refer to the data sheet or in the menu under "1-pt calibration".)

Sensor cartridge	Fresh air *	Standard gas: 21 vol.% O ₂ , 79 vol.% N ₂	N ₂	CO ₂
CO ₂ vol.% up to				
10	X	X		
25	X	X		
50		X	X	
50 biogas		X	X	
O ₂ vol.% up to				
25	X	X		
50	X	X		
100	X	X		
CH ₄ vol.% up to				
100			X	
100 biogas			X	
H ₂ vol.% up to				
10			X	X
50			X	X
100			X	X

*Fresh air means: 0.04 vol.% CO₂, 20.97 vol.% O₂, rest N₂

7.2. 2-point calibration

Notice

The 2-point calibration should only be carried out by experienced users! Please contact the **BlueSens gas sensor GmbH** service department for more information!

Notice

All previously performed 2-point calibrations of a cartridge can be reset in the sub-menu "Advanced settings".

Figure 68: First step in removing the cartridges

Figure 69: Last step in removing the cartridges

7.3. Filter change

The filters and sealing rings should be checked for damage and possible contamination (e.g. foam penetration) before each use and replaced if necessary. If there is continuous usage, we recommend changing the filter once a month.

7.3.1. Removing the sensor cartridges from the flow adapter

1. Turn the cartridge on the left side (you can read the label on cartridge 1) slightly to the right (anti-clockwise) and remove it from the flow adapter.
2. Repeat this until all sensor cartridges are removed.
3. Replace the filters and seals.

Notice!

Danger due to water entering the measuring cells of the **BlueVary**.

There is a risk of material damage and functional impairment.

→ The integrated filter is not intended for separating liquids (see section 4.2.1). Never install the sensor in such a way that liquid can enter.

If water enters the flow adapter, it must be blown out with a sufficient amount of dry air.

8. Appendix

8.1. Calibration table

The monthly 1-point calibration should be performed as described in section 7.1.

Fill in the following table during this procedure:


8.2. Pin assignment of the analog cable

art. no. Z-KA-00058

Colour	Pin assignment
White	4-20 mA output, RL < 500 Ohm for CH1
Brown	GND
Green	4-20 mA output, RL < 500 Ohm for CH2
Yellow	GND
Grey	4-20 mA output, RL < 500 Ohm for CH3
Pink	GND
Blue	1-point calibration, connect to GND for 5 s. *
Red	GND


* Only after it is activated in the Advanced Settings sub-menu, (see section 5.10).

8.3. Technical drawing

8.4. Technical data

8.4.1. BlueVary main unit

Electrical inputs and outputs	
Communication interfaces:	RS485 (Modbus RTU), Modbus TCP, RS232
Powersupply:	24 VDC
Maximum current consumption	1 A
Pin assignment of port A Snap-in 8-pin plug and M12 8-pin plug (Z-KA-00054) 	1. +24 V 2. GND 3. RS232_RXD 4. RS232_TXD 5. RS485_A+ 6. RS485_B- 7. RS485_GND 8. GND
The RS485 interface is galvanically isolated; therefore RS485_GND must be used.	

Electrical inputs and outputs	
Pin assignment Port B Snap-in 8-pin socket	1. CH1_out 2. GND 3. CH2_out 4. GND 5. CH3_out 6. GND 7. 1-point calibration, 5 s to GND* 8. GND
All 4-20 mA outputs are active outputs for RL < 500 Ohm. *Only after activation in the Advanced Settings sub-menu, (see section 5.10).	

Operating conditions	
Permissible ambient temperature range	5 – 55 °C during operations 0 – 60 °C during transport and storage
Permissible gas temperature	5 – 55 °C
Permissible ambient pressure range:	0.8 to 1.3 bar
Permissible relative pressure in exhaust hose	Max. 0.3 bar
Maximum duty cycle	100%

Mechanical specifications	
Width (incl. cartridges and ¼" flow adapter)	Approx. 130 mm
Height (incl. cartridges and ¼" flow adapter)	Approx. 265 mm (without cable and hose nozzles)
Depth (incl. cartridges and ¼" flow adapter)	Approx. 150 mm
Weight: BlueVary main unit only	Approx. 500 g
Weight: Complete system (BlueVary main unit, 3 cartridges and ¼" stainless steel flow adapter)	Approx. 2600 g
Protection class	IP65

Electrical input of power supply unit	
Voltage	100 - 240 V ~
Frequency	50 - 60 Hz
Maximum current consumption	0.7 A

8.4.2. **BlueVary** gas sensor cartridge

Refer to the separate gas sensor cartridge data sheet.

8.5. Warranty

The warranty period is regulated in the General Terms and Conditions.

Notice

The warranty claim may be voided if this operating manual is not observed, or in case of improper handling (e.g. opening the **BlueVary** of a sensor cartridge), or if a power supply unit other than the original one is used.

8.6. Service and support

Our qualified customer service is available to you as a partner.

If required, please contact your dealer or us directly:

Phone: +49 2366 4995 567 or e-mail to service@bluesens.de

USA: 847 201 3124 or e-mail to service@bluesens.com

8.7. Recycling and disposal

The **BlueVary** purchased from us, the accessories and the associated packaging are all recyclable.

- The **BlueVary** is packaged to protect it from transport damage. Please take any packaging that you no longer need to your local recycling centre.
- Recycle and dispose of the **BlueVary**, including its accessories and packaging, in accordance with the national legal requirements applicable.
- The  symbol on the **BlueVary** means that you are not allowed to dispose of the **BlueVary** with normal waste within the EU. It must be recycled separately.
- For disposal, we recommend returning the **BlueVary** and its accessories to **BlueSens gas sensor GmbH**. For this, you will need a **return authorisation number**, which you can obtain from our service team service@bluesens.de together with a valid return address. We guarantee the professional dismantling and environmentally friendly recycling of all components! This service is available free of charge for you.
- Any personal data on or in the **BlueVary** must be deleted before it is returned.

8.8. Legal notice

BlueVary operating manual - German

Index Z-BA-MULTI-BVARY

Version V2.00 Rev. 231220 002

© **BlueSens gas sensor GmbH**

The images and drawings in this operating manual may differ from the original; they are for illustrative purposes only. All information, including the technical information, is subject to change without notice. All images and graphics in this instruction manual: © **BlueSens gas sensor GmbH**, Herten.

The specifications in this manual have been carefully checked for accuracy. Nevertheless, **BlueSens gas sensor GmbH** assumes no liability for consequences arising from any errors in the description and illustrations. The General Terms and Conditions from **BlueSens gas sensor GmbH** apply.

BlueSens gas sensor GmbH

Snirgelskamp 25

D-45699 Herten, Germany

Phone: +049 (0)2366/4995500

Fax: +049 (0)2366/4995599

E-mail: service@bluesens.de

Home page: www.bluesens.de

Information status: December 2023

9. CE and FCC certificates

The **BlueVary** complies with the EMC Directive (2014/30/EU) using the harmonised standards EN 55011, EN 61326-1 and the RoHS Directive (2011/65/EU).

The Low Voltage Directive (2014/35/EU) does not apply because no voltage higher than 24 V is used.

9.1. CE CERTIFICATE

EU Declaration of Conformity according to

EU Directive 2014/30/EU (EMC Directive)

EU Directive 2011/65/EU (RoHS Directive).

complies with the essential requirements for the adoption of the laws of the Member States relating to electromagnetic compatibility set out in the Council Directive (2014/30/EU).

The assessment of the product for electromagnetic compatibility was carried out based on the following harmonised EC standards:

EMI EN 55011

EMS DE 61326-1

This declaration applies to all devices manufactured according to the tested prototype.

HERTEN, Germany 1 Mar. 2022

Location and date

Signature of manufacturer/
representative

We, the **BlueSens gas sensor GmbH**

Snirgelskamp 25, 45699 Herten (Germany)

hereby declare under our sole responsibility that the product:

Gas sensor

BlueVary including cartridges for CO₂, CH₄, O₂, O₂ec, H₂, pressure/humidity and pressure

when used together with our accessories:

Power supply,

Cable

9.2. FCC CERTIFICATE

Declaration of Conformity

We, the **BlueSens gas sensor GmbH**

Snriegelskamp 25, 45699 Herten (Germany)

hereby declare under our sole responsibility that the product:

This declaration applies to all devices manufactured according to the tested prototype.

HERTEN, Germany 1 Mar. 2022

Location and date

Signature of manufacturer/
representative

Gas sensor

BlueVary including cartridges for CO₂, CH₄, O₂, O₂ec, H₂, pressure/humidity and pressure.

when used together with our accessories:

Power supply,

Cable

for which this declaration is based, conforms to the following specifications:

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Wir, die **BlueSens gas sensor GmbH**