Operating manual

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YeastForce

MONITOR

BlueSens

.....

Content

1	ABOUT THIS DOCUMENT	3
1.1	Purpose	3
1.2	Target group	3
1.3	Symbols in use	3
2	FOR YOUR SAFETY	4
2.1	In general	4
2.2	Authorized personnel	4
2.3	Proper and intended usage	4
2.4	Warning against misuse	5
2.5	General safety instructions	5
2.6	CE conformity	5
3	PRODUCT DESCRIPTION	6
3.1	The sensor	6
3.2	Converting to a standardized gas	6
3.3	Accessories	7
3.4	Required laboratory materials	8
4	INSTALLATION	9
4.1	General instructions	9
4.2	Attaching the filter adapter	9
5	ELECTRICAL CONNECTION	10
5.1	In general	10
5.2	USB - RS232 cable wrap-around	10
5.3	RS232 connection to the PC	10
6	YEASTFORCE MONITOR	12
6.1	Installing the software	12
6.2	Connecting the sensors	14
6.3	Search for heads:	14
6.4	CO ₂ sensor: single-point calibration	16
6.5	Volume inspection	19
6.6	Preparing a measurement	22
6.7	Starting a measurement	26
7	APPENDIX	28
7.1	Technical data	28
7.2	Warranty	28
7.3	Service and Support	28
7.4	Proper disposal	29
7.5	Imprint	29
8	CE-CONFORMITY	30
9	FCC-CONFORMITY	31

1 About this document

1.1 Purpose

This operating manual provides you with all the information you need to quickly commission and safely operate the **YeastForce**. Please read this operating manual before you start the initial commissioning. Keep this operating manual in a safe place for future reference.

1.2 Target group

This operating manual is intended for trained specialist personnel. The contents of this manual must be provided to and implemented by the trained personnel.

1.3 Symbols in use

Danger!

This symbol indicates a possibly dangerous situation. Failure to comply with these safety instructions may result in personal injury.

Caution!

This symbol indicates a situation, which could result in damage to equipment or other property.

Note!

This symbol indicates additional helpful information.

1 Sequence of actions

Sequential numbering indicates successive steps.

2 For your safety

2.1 In general

The **YeastForce** was inspected before it left the factory and was in an operationally ready condition.

Please read this operating manual carefully before installing and commissioning the device. The operating manual contains safety instructions that must be followed to ensure safe operation.

The contents of this manual correspond to the status as of November 2020; they may be changed without prior notice. We reserve the right to make technical changes in the course of further developments.

The accuracy of the information in this manual has been carefully checked. Nevertheless, the BlueSens gas sensor GmbH assumes no liability for consequences arising from any errors in the description and illustrations. The General Terms and Conditions of BlueSens Gas Sensor GmbH apply.

This device must never be operated under conditions that do not conform to its specifications and information specified on the ratings plate.

Maintenance and repair may only be carried out by properly trained, competent personnel who are familiar with the associated risks and warranty provisions.

2.2 Authorized personnel

The operations described in this operating manual may only be carried out by trained specialist personnel who have been authorized by the facility operator. For safety and warranty reasons, any further intervention or repairs to this device may only be carried out by personnel from BlueSens gas sensor GmbH.

2.3 Proper and intended usage

The **YeastForce**, as described in the Technical Specifications, is a device for measuring volumetric flows in the specified flow range and under the specified conditions. It is used to monitor metabolic procedures within biological processes. The **YeastForce** measuring device may only be used in well-ventilated rooms.

**Danger!**

This measuring device has no ATEX approval and may only be used in well-ventilated rooms.

2.4 Warning against misuse

The YeastForce may not be used as a safety component for monitoring in facilities; it must not be used in explosive zones.

2.5 General safety instructions

This device may pose application-specific hazards if it is used improperly.

**Danger!**

Incorrect installation or configuration may result in an explosion.

Check all connections for leaks after the installation is finished.

2.6 CE conformity

The YeastForce complies with the EMC Directive (2014/30/EU) using the harmonized standards EN 55011, EN 61326-1.

The Low Voltage Directive (2014/35/EU) does not apply because no voltage higher than 24 V is used.

See the last pages of this manual for the CE and the FCC certificate.

3 Product description

3.1 The sensor

The **YeastForce** consists of a sensor top section [A] (with sensors for pressure, CO₂ concentration and temperature) and a lower section consisting of a jar [B] with one bracket [C].

There are two connectors on the top of the **YeastForce**. One connection is used for the power supply and communications [D]. The other connection is used for the gas discharge [E]. A filter is included for the gas discharge.

Figure 1

During operations, the filter must be mounted at the gas discharge outlet. Otherwise, the pressure cannot be released from the sensor. This would prevent the sensor from functioning correctly; dangerous over-pressure could then occur.

To achieve the best measurement result, the **YeastForce** should be operated in a bath of water at a constant temperature while measuring. The water level must not be higher than the glass jar.

3.2 Converting to a standardized gas

The **YeastForce** has sensors for pressure and temperature. With the help of these measured values, the internally calculated volume is normalized to the standard conditions of 1013.25 mbar and 273.15 K (0°C).

The calculation can be converted back to any laboratory conditions by using the general gas equation. Since the number of particles in a closed system is constant, the following formula applies:

$$\frac{p * V}{T} = \text{constant}$$

where

p = Pressure [mbar]

V = Volume [ml]

T = Temperature [K]

In the template settings of the YeastForce Monitor, you can set the standard conditions for the temperature to easily compare it with other measuring methods according to user requirements.

3.3 Accessories

The following accessories are required for the initial commissioning of the **YeastForce**:

- 3.3.1 **USB → RS232 serial port adapter with power supply connection**
(Z-KA-00014)

With wrap-around (refer to chapter 5.2)

Figure 2. Connection from the PC (USB) to the **YeastForce** (M12 / 8-pin)

- 3.3.2 **12 V DC power supply**
(Z-NT-00010)

Figure 3. 12 V DC power supply

- 3.3.3 **Glass jar for the YeastForce**
Standard empty volume of 562.5 ml (± 2 ml)
(Z-GF-00019)

Figure 4 Glass jar for the **YeastForce**

- 3.3.4 **Bracket for the YeastForce**
(Z-XX-00166)

Figure 5 Bracket

- 3.3.5 **Filter adapter for the YeastForce gas discharge**
(Z-XX-00133)

Figure 6 Filter adapter

- 3.3.6 **Piston pump for volume inspection (20 ml)**
(Z-XX-00134)

Figure 7 Piston pump

- 3.3.7 **Active USB hub**

Optional
(K-00006)

Figure 8 Active USB hub

- 3.3.8 **Seal ring**
(Z-XX-00200)

Figure 9 Seal Ring

3.4 Required laboratory materials

3.4.1 Sodium hydroxide NaOH

For removing residual CO₂ during the single-point calibration of the CO₂ sensor

Available from chemical specialty stores

3.4.2 Water bath

The sensor is calibrated to a temperature range of 15-40°C. Therefore the water bath may have a temperature range of 25-40°C.

Water height: 5 – 9 cm (maximum height of glass jar)

The **YeastForce** must not be completely immersed because it has a protection class of IP64!

4 Installation

4.1 General instructions

During delivery, the **YeastForce** has been shipped in protective packaging. This packaging protects against normal transport loads. Nevertheless, please check before installation whether the device has been damaged by improper transport or improper storage. Safe operation is not possible if the device is damaged. The device may not be installed or put into operation if damage is detected.

Make sure that the specified operating conditions are maintained at all times.

Installation should only be carried out under expert guidance and in accordance with the relevant recognized occupational safety regulations.

4.2 Attaching the filter adapter

Important!

If the sensor is supplied without a filter adapter on the gas discharge, the filter adapter must be attached to the gas discharge after unpacking the **YeastForce** sensors.

Figure 10 Attaching the filter adapter

5 Electrical connection

5.1 In general

Caution!

- Read the installation instructions carefully to avoid damaging this device.
- Proceed step by step.
- Use only the original plugs, cables and power supplies.
- Never connect or disconnect the plug when the device is connected to the power supply.
- This device has no on/off switch. It starts operating immediately after it is connected to the power supply.
- Incorrect operations may damage the device.

5.2 USB - RS232 cable wrap-around

To minimize the risk that the **YeastForce** cable tangles with any other cables, the power supply cable should be wrapped around the data cable. A simple method is to wrap the power supply cable around the data cable.

First, wrap the data cable around the M12 8-pin connector. Continue until the power supply cable is completely wrapped around the data cable.

At the end, make a small loop and form the cable into a knot. This prevents the power cable from unwinding.

Figure 11 Wrapping the cables

Figure 12 Knot the cables at the end

5.3 RS232 connection to the PC

The USB-to-RS232 cable (with power supply) is first attached to the **YeastForce** using the M12 8-pin adapter.

The plug of the power supply unit is then connected to the power connection of the data cable and secured so that it cannot slip out.

Figure 13

The power supply unit may now be plugged into a power outlet. Now, the sensor starts working.

The cable's USB plug should now be connected to a computer (in an active hub if necessary).

Figure 14

The sensor needs about one hour until it is ready for operations. No measurement may be made before this.

The **YeastForce Monitor** program shows a flashing icon during this warm-up period:

[Wait while sensor heats up]

Figure 15

6 YeastForce Monitor

You can download the latest software at the following link:

https://www.bluesens.com/fileadmin/user_upload/downloads/YeastForceMonitorSetup_1_4_20.msi

6.1 Installing the software

Double-click on the .msi file to start the installation. Since the program is installed in the Program Files folder, admin rights are required during installation.

The Welcome page appears.

- ➔ Click on [Next]


Figure 16

Select the installation path and user settings.

- ➔ Click on [Next]


Figure 17

Confirm and start the installation:

- ➔ Click on [Next]


Figure 18

The program is installed:

Please wait:

- A progress bar shows the installation progress.


Figure 19

The final dialogue box appears after a successful installation:

- Click on [Close]


Figure 20

6.2 Connecting the sensors

If the drivers are not found automatically, you can download the latest FTDI-VCP driver from the following page:

<https://www.ftdichip.com/Drivers/VCP.htm>

6.3 Search for heads:


Figure 21 **YeastForce Monitor**, after the program start

When the program is started, all available COM ports in the computer are scanned for connected sensors. These are then loaded into the program.

If not all sensors are found, the program can search for more sensors again via the [Heads→Search Heads] option.


Figure 22

If a found sensor is then disconnected and reconnected to the computer, it can be reactivated with the [Reconnect] option.


Figure 23

6.3.1 Configuring a sensor head

To configure the **YeastForce** sensor, the Configuration window can be invoked from the [Heads] menu or the [HeadInfo] right-click context menu.


Figure 24

Use the [Custom Head Color] drop-down list to change the display colour of the sensor.

The [Show Always] option causes the sensor to always be displayed on the user interface. If the sensor is not found, this is displayed immediately. You can reactivate the sensor in the program's context menu via [Reconnect].

If the [Show Export Options / Calibration] option is selected, additional information about the sensor is displayed.

This calibration information relates to the data about the operating status saved in the sensor or the program. If the single-point calibration or volume inspection was performed more than one month ago, it must be carried out again. The previous dates (of last execution) are displayed in this window.


Figure 25


Figure 26

6.4 CO₂ sensor: single-point calibration

CO₂ absorber NaOH

Danger!

NaOH may only be used in accordance with the applicable safety precautions!

Serious injury to skin or eyes

**Wear protective gloves
Wear safety goggles**

Gather the materials together:

- **YeastForce** sensor
- Glass jar for the **YeastForce**
- NaOH solution (e.g. 3 M ≈ 10 %)
- Small vessel for at least 20 ml NaOH (e.g. 100 ml beaker),
→ which will fit into the glass jar ←

Figure 27

Pour 20 to 100 ml NaOH into a suitable vessel and place it in the glass jar.

Figure 28

Attach the **YeastForce** to the jar using the bracket on the glass jar.

Figure 29

Start the single-point calibration in the **YeastForce Monitor** (via the right-click context menu for the sensor and confirm the information about the addition of NaOH by clicking on OK).


Figure 30

The concentration is then recorded for 5 minutes and the slope of the graph is displayed. The upper graph shows the progress over the entire time of the measurement. The lower graph always shows only the last five minutes. The slope of the graph is also calculated during these five minutes.


Figure 31

The CO₂ concentration now decreases until there are hardly any CO₂ molecules left in the vessel. If the slope drops below 2 ppm/min, the button for the single-point is enabled.

1-Point Calibration	
Slope:	-1.9 [ppm/min]
<input type="button" value="Perform 1-Point Calibration"/>	
<input type="button" value="Restart"/>	

Figure 32

After the single-point calibration is finished, the new concentration for the sensor is displayed. If this is more than 0.01% by volume, the single-point calibration should be performed again.


Figure 33

Click on the Restart button: the instrument waits again for 5 minutes until a new single-point calibration is possible.

Also here, the slope must be less than 2 ppm/min.


Figure 34

You can click on the [Return to YeastForce Monitor] button to return to the main program.


Figure 35

6.5 Volume inspection

Connect the sensor to the power supply. Connect the USB cable to the computer.

(Refer to chapter 5)

Mount the **YeastForce** on a glass jar. Place them together in a water bath with a water height between 5 and 9 cm.

The temperature in the water bath should be between 25 and 40 °C, and higher than room temperature.

Keep the water bath and sensor at a constant temperature for at least 30 minutes.

Figure 36

Start the **YeastForce Monitor** software program and search for the connected **YeastForce** head.

(Refer to chapter 6.3 Search for heads)

Start a volume inspection from the right-click context menu.


Figure 37

Pull out the piston pump.

Figure 38

Remove the **YeastForce**'s protective filter by pushing the fluted ring down towards the orange band.


Figure 39

Enter the temperature of the piston in the program. (here: 27.5°C)

If it is too cold, heat the piston to 20°C with your hands.

Figure 40

Place the piston pump on the **YeastForce** until the piston clicks into place.

Then, confirm this in the program:
[OK – Piston connected]


Figure 41

Wait about 10 seconds.

If the pressure is constant and the program prompts you to do so, press the piston rod into the piston slowly and consistently.


Figure 42

After the valve closes (audible soft click and the message [Remove Piston]), the piston can be removed and the filter can be put on.

Figure 43

The pressure rises in the **YeastForce**. The program records the pressure and calculates the tightness and volume of the system.

The lower graph shows a partial section of the upper graph.

This automatically displays the relevant part of the measurement with more details.

Figure 44

After about 5 minutes, the calculation is completed and the results are displayed in a window.


Figure 45

Click on the button [Show/Hide Details] to display the individual values of the calculation.

If the calculation was not successful, the inspection can be repeated by clicking the [Retry] button.


Figure 46

If the filter connector is not plugged back into the gas discharge, the sensor will report a message. The measurement should then be repeated to check the return for external pressure.


Figure 47

6.6 Preparing a measurement

6.6.1 Preparing a template


Figure 48

A template for a measurement can be configured for each type of dough.

The template configuration opens when you click on the context menu of a template field:


Figure 49

The volume of the glass jar and the volume of the dough are recorded with an accuracy of ± 2 ml when possible.

The total volume (Total [ml]) of the system >>YeastForce Sensor + jar << and the dead volume are calculated based on the entered values and the additional sensor volume (~ 66 ml).

Four characteristics of the dough must be specified in order to categorize the different types of dough.

Refer to chapter 6.6.2.


Figure 50

The selection options can be changed, created or deleted by the user with the [*<New>* Or *<Edit>*] menu option.


Figure 51

The measurement duration, the logging interval and the calculation temperature can be set in Measurement Options. A start delay is currently deactivated. The calculation temperature allows the standard temperature to be set to a different temperature so that a better comparison can be made with other measuring methods.

The [Log Raw Data] option may be enabled so that the data can be subsequently evaluated by **Blue-Sens**. This means that all information received by the program every second is saved in a separate file.

You can create separate settings for the report for each template.

The [Summary Interval] indicates the time interval for the measured values in the summary. The display of the graphs can be changed so that the dough volume and the gas holding capacity are not shown.

The [Create PDF] option causes the report to be saved as a PDF file in the Reports folder.

The key information for each measurement is normally stored in a report log in .csv format. If the corresponding option is not selected, no entry is created.

The user may specify a default comment. There are two options for this:

The [Use Default Comment] option ensures that the specified comment is added to the sensor information when the template is selected.

If the upper option is selected, the measurement cannot be started until the user has changed the comment (e.g. to enter additional information such as a batch number).


Figure 52


Figure 53


Figure 54

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The user may edit each entry and add user-defined entries.

6.6.3 Creating a new template

There are three different methods to create a template in the program:

Via the right-click context menu for a template, an existing template can be copied to a new template [Copy Into New Template]. You only need to specify a new name. Then, the template can be changed to fit your requirements.


Figure 55

Click on the menu option [Templates→New] to create a new template with standard settings.

A template can be loaded from a template*.xml file by selecting the menu option [Templates→Load]. Information on exchanging templates can be obtained from **BlueSens**.


Figure 56

6.7 Starting a measurement

Before starting a measurement, a template must first be selected in the [Selection] tab. A selected template is then framed in the colour of the sensor.

Figure 57

As soon as a template field is selected, the cursor jumps to the comment field. The comment may be changed but must not be empty.


Figure 58

If error messages are present, no new measurement can be started. If an attempt is made to start a measurement nevertheless, an error message is displayed.


Figure 59

When there is no error message, the measurement can be started by clicking on the Start symbol.

The Start symbol is then disabled and the Stop symbol is enabled.

Figure 60

After a measurement is started, the [Chart] tab is enabled and the graph of the volume produced is displayed.

The display for the [dough volume] and [gas retention] can be selected from tabs at the bottom.

6.7.1 Status icons

Icon	Icon 2 (for GIFs)	[Error Text] Information
		<p>[Wait while sensor heats up] Wait until the sensor has heated up.</p> <p>→ This is displayed for approx. 1 hour after connecting the sensor to the power supply.</p>
		<p>[Signal too low for measuring] The sensor signal is too low.</p> <p>If this error does not disappear within 20 minutes, the sensor lamp may be defective.</p> <p>→ Please contact the service department.</p>
		<p>[The Last Volume Inspection has not been successful] The previous volume inspection was not successful.</p> <p>→ Please repeat the volume inspection. (refer to 6.5)</p>
		<p>[The last Volume Inspection dates back more than one month] The last volume inspection was made more than a month ago.</p> <p>→ Please carry out a volume inspection: (refer to 6.5)</p>
		<p>[The last 1Point Calibration dates back more than one month] The last single-point calibration was made more than a month ago.</p> <p>→ Please carry out a single-point calibration: (refer to 6.4)</p>
		<p>[The Connection to the Sensor is Lost] The connection to the sensor is no longer established.</p> <p>→ Check the data cable and power supply.</p>
	(rotating)	<p>The concentration is too high to start the measurement. Please flush the sensor out with fresh air.</p> <p>(This appears as long as the concentration in the sensor at the start of the measurement is > 0.5 vol.-% CO₂)</p> <p>→ Rinse sensor with fresh air = Place sensor on its side.</p> <p>→ If that is inadequate: Single-point calibration?</p>
		<p>Open the report PDF file of the previous measurement</p> <p>Open the report log file in csv format</p>

7 Appendix

7.1 Technical data

Number of measuring channels	3: CO ₂ , pressure, temperature
Measuring method for CO ₂	Two infra-red wavelengths
Pressure range	0.8 – 1.3 bar
Accuracy	< 0.2 % full scale ± 3 % reading
Measuring range for CO ₂	25 % by volume
Measurable gas volume flow	0 – 360 ml/h in 500 ml sample vessel
Resolution for volume and CO ₂	V = 0.01 ml; C = 0.005 vol.%
Dimensions	Ø 115 x 205 mm
Temperature range	15 – 40 °C
Ambient humidity	< 90 %
Housing	IP64
PC connection	USB
Measuring interval	1 s (internal), Logging interval (default) 10 s, Summary interval (default) 10 min
Electrical supply	12 V DC via power supply unit (included)
Requirements for the PC	8 GB of RAM, Operating system: Win7 64 bit or later, One USB2 port, Display resolution > 1280 x 800 px
Software	YeastForce Monitor
Sample weight	Depending on the size of the sample vessel, Volume ratio 1:10
Size of the sample vessel	10 times larger than dough sample/volume

7.2 Warranty

The warranty period for **YeastForce** is stated in our terms and conditions.

A loss of the warranty is likely to be caused by non-observance of these operating instructions, in case of improper handling (e.g. opening the **YeastForce**) and when using an other than the supplied original power supply.

7.3 Service and Support

Our qualified customer service will be happy to assist you as a partner.

If necessary, please contact your distributor or directly to us:

Tel.: +49 2366 4995 567 or via e-mail service@bluesens.de

USA: 847 201 3124 or via e-mail service@bluesens.com

7.4 Proper disposal

We recommend sending back old or broken **YeastForce** to BlueSens gas sensor GmbH (address see under "imprint"). We will then guarantee you a professional decomposition and environmentally compatible recycling of all components!

7.5 Imprint

YeastForce Operating Manual English

Edition Rev.220301 003

© BlueSens gas sensor GmbH

The pictures and drawings in this manual may differ from the original, they are only for illustration purposes. All specifications – including technical specifications – can be changed without notice. All pictures and graphics in this guide: © BlueSens gas sensor GmbH, Herten

BlueSens gas sensor GmbH

Snirgelskamp 25

D-45699 Herten, Germany

Phone +049 (0)2366/4995500

Fax +049 (0)2366/4995599

E-mail: service@bluesens.de

- State of information: March 2022 –

8 CE-CONFORMITY

EU Declaration of Conformity according
EU Directive 2014/30/EU (EMC - Directive)
EU Directive 2011/65/EU (RoHS-Directive)

We, BlueSens gas sensor GmbH
Snirgelskamp 25, 45699 Herten (Germany)
herewith declare in sole responsibility that the product:

Sensor
YeastForce

In combination with our accessories:

Power supply,
Cables

observe the essential requirements determined in the Council Directive for
the adoption of the legal regulations of the Member States about electro-
magnetic compatibility (2014/30/EU).

The judgement of the product to electromagnetic compatibility was effected
on basis of the following EC harmonised standards:

EMI EN 55011
EMS EN 61326-1

This declaration applies to all specimen manufactured according to the
sample tested.

HERTEN, 20.10.2019

Place and Date
representative

Countersign of manufacturer/rep-

9 FCC-CONFORMITY

Declaration of Conformity

We, BlueSens gas sensor GmbH

Snirgelskamp 25, 45699 Herten (Germany)

herewith declare in sole responsibility that the product:

Sensor

YeastForce

In combination with our accessories:

Power supply,

Cables

To which this declaration relates is in conformity with the following specifications:

FCC 47 CFR, 15B, §15.107 & §15.109

This declaration applies to all specimen manufactured according to the sample tested.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

HERTEN, 20.10.2019

Place and Date

Countersign of manufacturer/representative

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Wir, die BlueSens gas sensor GmbH