Updating Manual multilingual

BlueVary

BlueSens
.....

Updating Manual

BlueVary

BlueSens
.....

Content

| | | |
|-------|--|----|
| 1 | ABOUT THIS DOCUMENT | 3 |
| 1.1 | Function..... | 3 |
| 1.2 | Target group..... | 3 |
| 1.3 | Symbols used..... | 3 |
| 2 | FOR YOUR SAFETY..... | 4 |
| 2.1 | General information..... | 4 |
| 2.2 | Authorized personnel | 4 |
| 2.3 | Proper use..... | 4 |
| 2.4 | Misuse warning | 4 |
| 2.5 | General safety instructions..... | 4 |
| 3 | INSTALLATION..... | 5 |
| 3.1 | General instructions | 5 |
| 3.2 | Electrical connection | 5 |
| 3.2.1 | General information..... | 5 |
| 3.2.2 | Power up | 6 |
| 3.2.3 | Connection via RS232/USB..... | 6 |
| 4 | PREPARATION OF THE BLUEVARY VIA THE MENU | 7 |
| 4.1 | General explanations | 7 |
| 4.2 | Starting display and main menu..... | 8 |
| 4.3 | Selecting Admin user | 9 |
| 4.4 | Switch BlueVary to update mode | 10 |
| 5 | UPDATING THE FIRMWARE | 11 |
| 5.1 | Preparing the software | 11 |
| 5.2 | Performing the update..... | 11 |
| 6 | APPENDIX | 15 |
| 6.1 | Warranty..... | 15 |
| 6.2 | Service and Support..... | 15 |
| 6.3 | Imprint | 15 |

1 About this document

1.1 Function

This manual provides you with all of the necessary information for updating the firmware of the **BlueVary**. Therefore, please read this manual before starting a firmware update. Keep this manual safe for later usage.

1.2 Target group

This manual is intended for use by trained specialist personnel. The content of this manual must be made available to the corresponding personnel and followed by them.

1.3 Symbols used

Hazard!

This symbol indicates a possibly hazardous situation. Failure to observe this safety instruction can result in personal injury.

Caution!

This symbol indicates the possibility of damage to property.

Note!

This symbol indicates helpful, additional information.

List

The bullet point indicates a list in which the order of items has no particular significance.

Sequence of actions

Numbers indicate steps to be performed in a certain sequence.

2 For your safety

2.1 General information

The **BlueVary** was programmed with the latest firmware version when it left our plant.

Please read this manual carefully before updating the firmware. The manual contains safety instructions that must be observed to ensure safe operation.

The device may never be operated in conditions that do not correspond to the specifications on the specification plate.

Maintenance and servicing may only be performed by specially trained personnel who are familiar with the hazards inherent to the work as well as the guarantee terms.

2.2 Authorized personnel

All actions described in this manual may only be conducted by specialist personnel who have been authorized by the plant operator. Work on the device other than that described in this manual and the operation manual may only be performed by personnel of the BlueSens gas sensor GmbH company for safety reasons and to ensure compliance with the terms of the guarantee.

2.3 Proper use

The **BlueVary** is a gas sensor for measuring gas concentrations for the specified gases in the specified concentration area under the conditions described in the technical data. It is used to monitor metabolism in biological processes such as fermentation.

2.4 Misuse warning

The **BlueVary** may not be used as a safety component for monitoring gases in systems or be used as a gas warning device.

2.5 General safety instructions

If the device is not used correctly, or used for its intended purpose, application-specific hazards may arise.

Danger!

**If the device is incorrectly installed or set, there is a danger of explosions and poisoning.
After installation, check all connections for leaks.**

Note! Please note the safety instructions in all the following chapters, especially in the chapters 3 and 4 of the user manual.

3 Installation

3.1 General instructions

The **BlueVary** was protected on its way to the installation location by a special packaging (Keep this package for maintenance purposes). This secures it against the usual transport strains. However, before installation, check if the device has suffered damaged due to improper transport or incorrect storage. If the device is damaged in any way, operation without hazards is not possible and the device may not be installed and taken into operation.

Check whether the enclosed materials, such as gaskets, are suitable for your process conditions (pressure, temperature, etc.).

The installation should only be performed under supervision by a specialist and in compliance with all applicable work safety rules.

3.2 Electrical connection

3.2.1 General information

Caution!

Read the installation instructions carefully to avoid damage to the device.

Proceed step-by-step.

Only use the original plugs, cables and power adapters (Fig. 1).

Never connect or disconnect plugs or cartridges when the device is connected to the power supply.

The device does not have an on/off switch; it starts operation as soon as it is connected to the power supply.

Improper operation can result in damage to the device.

Fig. 1: original cables and electrical accessories.

3.2.2 Power up

1. Make your choice of the corresponding blade B (fig. 1) necessary in your country to plug in the supplied power supply A (fig. 1).
2. Connect cable E (fig.1) to port A (fig. 2).
3. Connect cable E (fig. 1) and cable C (fig. 1) (use connector with the two wires)
4. Attach power supply to cable C (fig.3)

Fig. 2: Port A: Power in, RS232, ModBus RTU
Port B: 3 x Analogue outputs


Fig. 3: Connection of power supply A to connector of cable C

3.2.3 Connection via RS232/USB

The following items will be needed for the electrical connection RS232 to USB (fig. 4):

1. Connect cable D to the corresponding connector of cable C
2. Plug in the USB adapter into a free USB port of your computer. Be aware that the USB adapter does not replace the power supply!
3. After switching on the computer a COM port should be automatically created in Windows. Some Windows versions require a driver. Check the USB flash drive (if BlueVis is purchased) or visit the FTDI website for the latest RS232 driver.


Fig. 4: RS232-sticker on USB adapter of cable D

4 Preparation of the BlueVary via the menu

4.1 General explanations

The **BlueVary** has an integrated keypad and a LCD display for showing values and changing settings. To indicate the selected line a „-“ [cursor] is displayed in front of the text of the line (see fig. 5).

The functions of the keypad are:

↑: Moves the curser upwards / increases values

↓: Moves the curser downwards / decreases values

←: Moves the curser to the left

→: Moves the curser to the right

OK: Input acknowledgement / selections

esc: End of input / back to last level


Fig. 5: Available buttons on the **BlueVary** Center. ↑↓ ←→ OK esc

Glossary of the pictures: **Bold** type lines are for input or to make a choice.

The keypad is **Bold** type within the text: ↑↓ ←→ OK esc.

Menu items are displayed *italic* within the text.

To select a menu item move the cursor with ↑↓ to the corresponding line and press **Ok**

For options, these are selected via ↑↓, numbers are entered via ↑↓ ←→

To finish an input press the **esc** – button.

To store values, it is necessary to select the menu item Apply settings.

Should a menu to be left, always can be pressed the **esc** – button.

4.2 Starting display and main menu

The starting display is shown in Fig. 6. It shows the gas type and concentration of the attached cartridges followed by the pressure and the absolute humidity. It also shows the lifetime of the attached cartridges.

Fig. 6: Starting display, schematic and picture

Pressing **OK** will lead you to the main menu:


Fig. 7: Starting display and the items of the main menu (schematic)

4.3 Selecting Admin user

In the submenu *Set user level* the access rights can be chosen according to the following list:

| user level | Access rights | Code |
|--------------------------|--|-------------|
| <i>Normal user</i> | Reading rights | Unnecessary |
| <i>Calibration user</i> | Normal user + calibration rights (1-point / 2-point) | 0001 |
| <i>Installation user</i> | Normal user + configuration rights | 2020 |
| <i>Admin</i> | Unlimited rights | 1818 |

☞ **After a reset or power up user level always switches to *Normal User*!**

To prepare the **BlueVary** for the firmware update the *Admin* has to be selected. Place the cursor with the $\uparrow\downarrow$ buttons in front of the *Admin*. By pressing **OK** button the user is selected. Now proceed with code input. Select the *Code* line ($\uparrow\downarrow$ + **OK**) and insert the code 1818 ($\uparrow\downarrow\leftarrow\rightarrow$ + **ESC**).


Fig. 8: Schematic of the submenu user level

To accept the input go to *Apply settings* ($\uparrow\downarrow$) and press **OK**.

Depending on whether the code was right or not one of the following messages appear:

Fig. 9: Possible messages after trying to apply settings in the User menu

The message has to be confirmed by pressing **OK**.

4.4 Switch BlueVary to update mode

The submenu *Advanced settings* gives the user access to some additional functions of the **BlueVary** and is shown in Fig. 10.


Fig. 10: Schematic of the submenu *Advanced settings*

To finish the preparation of the **BlueVary** for the firmware update the item *Read new firmware* has to be selected ($\uparrow\downarrow$) and press **OK**.

On the new screen that comes up (Fig. 11) **Reset** has to be selected ($\uparrow\downarrow$) and press **OK**.


Fig. 11: Schematic of the screen after selecting „Read new firmware“

The **BlueVary** now will execute a reset to start the firmware update mode.

☞ Be aware, that display may freeze or switch off. That is normal.

☞ The menu settings of the **BlueVary** are retained after a firmware update.

5 Updating the firmware

Take care that the **BlueVary** is connected to the power supply and via the RS232-USB-adaptor to the computer (Chapter 3). If the **BlueVary** is connected to a program via the adaptor, it has to be closed. When the communication works with a background service, also the service has to be closed. Always use the latest firmware and the latest version of the update program BlueSensBoot.exe. Both can be received from the BlueSens service team.

5.1 Preparing the software

The firmware and the update program are normally provided as a zip-file containing a folder with the latest firmware and the update program **BlueSensBoot.exe**. Unzip the File **BlueVary FirmwareUpdater.zip** to your desktop.


Fig. 12: Folder with Software (highlighted) and firmware

5.2 Performing the update

To start the software double-click on **BlueSensBoot.exe** (Fig.12).

After the software has started click on the *choose program* button (Fig 13) and select the BlueVary_VX_Y*.bin firmware file which is in the same folder as the program.


Fig. 13: Update program after starting

The next step is to select the right COM-Port. If the port number of the RS232-USB adaptor of the **BlueVary** is unknown, all available port numbers that are shown in the port list of BlueSensBoot (Fig. 14) must be written down. Then the program has to be closed and the adaptor disconnected from the computer. After opening the program again, the now missing port number has to be marked. After closing the program again, reconnecting the adaptor to the same USB port as before and opening the program again this port has to be chosen.


Fig. 14: List of available COM-Ports

If the **BlueVary** is in firmware update mode (see Chapter 4), it now can be connected to the program by clicking the *connect to device* button. (Fig. 15)


Fig. 15: Click on the connect to device button

When the **BlueVary** has been switched to the firmware update mode and the right COM port has been selected, a message will appear in the software and the programming can be started by clicking the *start programming* button. (Fig. 16)


Fig. 16: When message states the connection start programming can be clicked

Caution!

During the programming procedure the **BlueVary** must not be disconnected from computer or power supply. This can destroy the **BlueVary**.

The updating of the firmware will take less than 2 minutes. It is performed in two different steps. First the old firmware will be erased. (Fig. 17)


Fig. 17: First step, old firmware is erased

Then the new firmware is uploaded. (Fig. 18)


Fig. 18: Second step, new firmware is uploaded

When the programming is completed the program will switch the **BlueVary** from bootloader (firmware update) to application mode and restart it. To free the COM port for other application the **Disconnect** button has to be pressed. (Fig. 19)


Fig. 19: Software switched **BlueVary** to application mode, **Disconnect** button has to be pressed

Now the program has to be closed or another **BlueVary** has to be connected to the RS232-USB adaptor.

Please be aware, that the **BlueVary might have to heat up again.
(Yellow backlight of the display is on)**

6 Appendix

6.1 Warranty

The warranty period for **BlueVary** is stated in our terms and conditions.

A loss of the warranty is likely to be caused by non-observance of these instructions and the operating instructions, in case of improper handling (e.g. opening the **BlueVary** or the sensor cartridges) and when using an other than the supplied original power supply.

6.2 Service and Support

Our qualified customer service will be happy to assist you as a partner.
If necessary, please contact your distributor or directly to us:

Tel.: +49 2366 4995 567 or via e-mail service@bluesens.de

USA: 847 201 3124 or via e-mail service@bluesens.com

6.3 Imprint

BlueVary Update Manual English

Edition Rev 200421 001

© BlueSens gas sensor GmbH

The pictures and drawings in this manual may differ from the original, they are only for illustration purposes. All specifications – including technical specifications – can be changed without notice. All pictures and graphics in this guide: © BlueSens gas sensor GmbH, Herten

BlueSens gas sensor GmbH
Snirgelskamp 25
D-45699 Herten, Germany
Phone +049 (0)2366/4995500
Fax +049 (0)2366/4995599
E-mail: service@bluesens.de

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